


Connector type




Sensor type


S0013A Series Adult Disposable SpO₂ Sensors
Application site : Adult > 30kg,
Application site: Index or other finger




S0013A Series Pediatric Disposable SpO₂ Sensor
Application: Pediatric >10-50kg
Application site: Index or other finger



S0013A Series Infant Disposable SpO₂ Sensors
Application: Infant 3~20kg
Application site: great toe



S0013A Series Neonatal Disposable SpO₂ Sensors
Application: Neonatal <3 kg,
Adult >40kg,
Application site: Neonatal foot,
Adult Index or other finger




Medical grade application tapes



Microfoam (standard)



Transpore



Medaplast





Notice: for Transpore or Medaplast tape, please indicated on your orders as follows: T=Transpore, e.g. S0010J-ST; P=Medaplast, e.g.S0010J-SP

Type Disignation

S□□□□□-□

- Cable Length: L: indicated L=0.9m, S:indicated L=0.45m
- Disposable Probe type: J: indicated Adult disposable SpO₂ sensor; K: pediatric; L: infant; M: neonatal
- Brand code: 0010/Nellcor, 0036/Datex
- S: Indicated SpO₂ sensor

Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	OEM P/N	MED-LINK REF NO.	SENSOR TYPE	EMITTER	PRICE CODE
BCI (Biochem International)	BCI(Biochem International): 3401,3304,3303,3302,3301,3300; Spacelabs Burdick: PRO2; Dräger: Dräger Oxipac 2000	050 		S0076J-S	Adult Disposable SpO2 Sensor	660/905-2	A5
				S0076K-S	Pediatric Disposable SpO2 Sensor	660/905-2	A6
				S0076L-L	Infant Disposable SpO2 Sensor	660/905-2	A6
				S0076M-L	Neonatal Disposable SpO2 Sensor	660/905-2	A5
CSI Criticare Systems International	Alaris(IVAC):VitalCheck 4400,CSI:503DX,503Spot,504DX,506, 506DX,507ELC2 ScholarIII,8100Poet Plus,Poet (602-3A) (Sensor + Cable),Scholar (Sensor + Cable)	064 		S0044J-S	Adult Disposable SpO2 Sensor	660/940-3	A5
				S0044K-S	Pediatric Disposable SpO2 Sensor	660/940-3	A6
				S0044L-L	Infant Disposable SpO2 Sensor	660/940-3	A6
				S0044M-L	Neonatal Disposable SpO2 Sensor	660/940-3	A5
DATASCOPE	Datascope:Accutorr 3/4 SAT, Accutorr Plus,Passport (Database SpO2),Compatible for 0020-00-071-01	129 		S0052J-S	Adult Disposable SpO2 Sensor	660/940-3	A5
				S0052K-S	Pediatric Disposable SpO2 Sensor	660/940-3	A6
				S0052L-L	Infant Disposable SpO2 Sensor	660/940-3	A6
				S0052M-L	Neonatal Disposable SpO2 Sensor	660/940-3	A5
DATEX 德恩	AS/3,CS/3,CH/S,CH/2S,CH-RS,CG/2S,Cardiocap I,Cardiocap II,Cardiocap 5,Capnomac Ultima,Light Monitor,M-Estpr,Oscar II,OSP-200,Satlite Trans,Satlite,Trans Plus,System 5,ULTS00,ULTSIO,ULTSVI,ULTSVO, original factory code :SAF-F4	058 		S0036J-S	Adult Disposable SpO2 Sensor	660/905-2	A5
				S0036K-S	Pediatric Disposable SpO2 Sensor	660/905-2	A6
				S0036L-L	Infant Disposable SpO2 Sensor	660/905-2	A6
				S0036M-L	Neonatal Disposable SpO2 Sensor	660/905-2	A5
OHMEDA 德恩欧美达	Datex-Ohmeda:3740,3710,4700 OxiCap,5250RGM, Biox3760,Ohmeda3700,3770,3775,Rascal II, GE:7000Series,Corometrics511, TRAM A.R, TRAM Modular x00SL Series, as200SL,400SL,600SL,800A,800SL,Mennen :Horizon XL800 (Ohmeda SpO2)	057 		S0014J-S	Adult Disposable SpO2 Sensor	660/940-3	A5
				S0014K-S	Pediatric Disposable SpO2 Sensor	660/940-3	A6
				S0014L-L	Infant Disposable SpO2 Sensor	660/940-3	A6
				S0014M-L	Neonatal Disposable SpO2 Sensor	660/940-3	A5
OHMEDA 欧美达	Datex-Ohmeda:3740, 3710, only use for Ohmeda oximeter, EMITTER: 660/940-3			S0123J-S	Adult Disposable SpO2 Sensor	660/940-3	A9
				S0123K-S	Pediatric Disposable SpO2 Sensor	660/940-3	A0
				S0123L-L	Infant Disposable SpO2 Sensor	660/940-3	A0
				S0123M-L	Neonatal Disposable SpO2 Sensor	660/940-3	A9
DOLPHIN (OSI) 海豚	Dolphin ONE Oximeters, *digital PD output directly DS. Copatible original : 360	131 		S0099J-S	Adult Disposable SpO2 Sensor	660/905-2	A5
				S0099K-S	Pediatric Disposable SpO2 Sensor	660/905-2	A6
				S0099L-L	Infant Disposable SpO2 Sensor	660/905-2	A6
				S0099M-L	Neonatal Disposable SpO2 Sensor	660/905-2	A5
Masimo (original) 马斯默 (原装)	Redical	115 		S0020J-S	Adult Disposable SpO2 Sensor	660/905-2	D2
				S0020K-S	Pediatric Disposable SpO2 Sensor	660/905-2	D2
				S0020L-L	Infant Disposable SpO2 Sensor	660/905-2	D2
				S0020M-L	Neonatal Disposable SpO2 Sensor	660/905-2	D2
NELLCOR 奈普特	Nellcor®:N100 (req. preamp cable),N180,N200 (req. preamp cable),N-20,NPB-40,NPB-75,NPB-290/5 ,N-390,N-395,N-3000,N-6000	059 		S0010J-S	Adult Disposable SpO2 Sensor	660/905-2	A5
				S0010K-S	Pediatric Disposable SpO2 Sensor	660/905-2	A6
				S0010L-L	Infant Disposable SpO2 Sensor	660/905-2	A6
				S0010M-L	Neonatal Disposable SpO2 Sensor	660/905-2	A5
NELLCOR OXIMAX TECH (ORIGINAL) 奈普特 (原装)	All Nellcor & Nellcor equipped OEM monitors not listed above	049 		S0026J-S	Adult Disposable SpO2 Sensor	660/905-2	C2
				S0026K-S	Pediatric Disposable SpO2 Sensor	660/905-2	C2
				S0026L-L	Infant Disposable SpO2 Sensor	660/905-2	C2
				S0026M-L	Neonatal Disposable SpO2 Sensor	660/905-2	C2
NIHON KOHDEN 日本光电	AL-800PA Module,BSM-1101 LifeScope EC, BSM-2101/2 LifeScope L, BSM-3101 LifeScope LC, BSM-7100/7200 LifeScope8, BSM-8200/1 LifeScope 12, BSM-8301/2 LifeScope 9, SM-8800 LifeScope14,OLV-1100A	057 		S0069J-S	Adult Disposable SpO2 Sensor	660/940-3	A5
				S0069K-S	Pediatric Disposable SpO2 Sensor	660/940-3	A6
				S0069L-L	Infant Disposable SpO2 Sensor	660/940-3	A6
				S0069M-L	Neonatal Disposable SpO2 Sensor	660/940-3	A5

Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	OEM P/N	MED-LINK REF NO.	SENSOR TYPE	EMITTER	PRICE CODE
NIHON KOHDEN 日本光电	Bedside Monitor: BSM-4101,4103, 4111,4113 Life scope P, BSM-5105,5135 Life scope A, BSM-2351,2353 Life scope L, BSM-2301,2303 Life scope I			S0005J-S	Adult Disposable SpO ₂ Sensor	660/940-3	A5
				S0005K-S	Pediatric Disposable SpO ₂ Sensor	660/940-3	A6
				S0005L-L	Infant Disposable SpO ₂ Sensor	660/940-3	A6
				S0005M-L	Neonatal Disposable SpO ₂ Sensor	660/940-3	A5
NONIN 燕牌	NONIN: 8500,8600,8700,8800			S0080J-S	Adult Disposable SpO ₂ Sensor	660/905-2	A5
				S0080K-S	Pediatric Disposable SpO ₂ Sensor	660/905-2	A6
				S0080L-L	Infant Disposable SpO ₂ Sensor	660/905-2	A6
				S0080M-L	Neonatal Disposable SpO ₂ Sensor	660/905-2	A5
NOVANMETRIX 诺医	Novamatrix:505,510,511,515B/C,5207000,7100,Oxypleth; original match : 8776-00; Spacelabs:90351D/G,90465,90466,90467,90489,90651-07/08,90419-01,90651A-07,original match:8776-00			S0062J-S	Adult Disposable SpO ₂ Sensor	660/940-3	A5
				S0062K-S	Pediatric Disposable SpO ₂ Sensor	660/940-3	A6
				S0062L-L	Infant Disposable SpO ₂ Sensor	660/940-3	A6
				S0062M-L	Neonatal Disposable SpO ₂ Sensor	660/940-3	A5
PALCO 派克	Palco: ALL			S0079J-S	Adult Disposable SpO ₂ Sensor	660/905-2	A5
				S0079K-S	Pediatric Disposable SpO ₂ Sensor	660/905-2	A6
				S0079L-L	Infant Disposable SpO ₂ Sensor	660/905-2	A6
				S0079M-L	Neonatal Disposable SpO ₂ Sensor	660/905-2	A5

SPECIFICATIONS

Oxygen Saturation Range (%SpO₂) : 0% to 100%

Pulse Rate Range(bpm) : 30 to 245

Accuracy

%SpO₂) (±1 S.D.)

90-100%±2 digits, 70-89% ± 3 digits, 0-69% unspecified)

Pulse Rate: ±3 bpm

LEDS:

Red: 660nm nominal,

Infrared: 880/905/940nm nominal

Storage and Transport Environment: -10 °C~+40 °C

Operating Temperature: +5 °C~+40 °C

Humidity

Operating/Storage/Transportation 0% to 80% non condensing Altitude

Hyperbaric Pressure 86kPa~106kPa

Biocompatibility:

All patient contact materials are latex free and have been tested according to ISO 10993-5,ISO 10993-10.

MED-LINKET is pleased to introduce multiple sensor application tapes to better meet your needs, the sensor product line gives you the flexibility and the choice in applying sensors utilizing medical grade medaplast application tapes.

The MED-LINKET line of disposable sensors feature a unique pre-aligned design that ensures accurate sensor placement, as well as a means to ensure sensor functionality. The application tapes maximize patient comfort and allow easy sensor repositioning.

Other compatible SpO₂ sensors are available and are compatible with most commonly used oximeter products and modules. MED-LINKET complete line of sensors are designed to provide high quality performance and are affordably priced.