









Disposable ECG Electrodes





Compatible Brand	Product Guide	OEM P/N	MED-LINK REF NO.	Description	PRICE CODE
all Monitoring/ Telemetry/ Electrocardiograph instrument		GEMS OEM P/N: E9001BG	V0014C/P	Disposable Paediatric ECG Electrodes, Cloth Round, 40mm diameter, wet gel, snap electrode. Silver/silver chloride (Ag/AgCl) sensor, repositionable. - packaging (film/disc): 1000	J0
		GEMS P/N: E9007SG	V0014A	Disposable Adult ECG Electrodes, Cloth Round, 2" (50mm) diameter, wet gel, snap electrode. Silver/silver chloride (Ag/AgCl) sensor, repositionable. - packaging (film/disc): 1000	J0

Cardiography Reuseble ECG Electrodes 心电图仪可重复使用心电电极



● Suction ECG electrodes 胸吸球心电电极

Compatible Brand	Product Guide	OEM P/N	MED-LINK REF NO.	Description	PRICE CODE
all Resting ECG Recommended for Nihon Konden, Fukuda, NEC, Kenz, China		Nihon konden: H043	V0013A	Adult suction electrode, - perfect fit for <u>3.0mm Din plug</u> leadwires - Latex-free, - Adult size replacement bulb for all 21mm diameter electrodes. Bulb diameter is 26mm. - Metal Electrodes: Nickel-plate. - 1 box = 6 bulbs	A0
			V0013B	Adult suction electrode, - perfect fit for <u>0.156-in.(4mm) banana plug</u> leadwires - Latex-free - Adult size replacement bulb for all 21mm diameter electrodes. Bulb diameter is 26mm. - Metal Electrodes: Nickel-plate. - 1 box = 6 bulbs	A0
all Resting ECG Recommended for Nihon Konden, Fukuda, NEC, Kenz, China			V0013A-P	Paediatric suction electrode, - perfect fit for <u>3.0mm Din plug</u> leadwires - Latex-free, - Adult size replacement bulb for all 17mm diameter electrodes. Bulb diameter is 23mm. - Metal Electrodes: Nickel-plate. - 1 box = 6 bulbs	A0
			V0013B-P	Paediatric suction electrode, - perfect fit for <u>0.156-in.(4mm) banana plug</u> leadwires - Latex-free - Adult size replacement bulb for all 17mm diameter electrodes. Bulb diameter is 23mm. - Metal Electrodes: Nickel-plate. - 1 box = 6 bulbs	A0
all Resting ECG Recommended for GE-Marquette, Schiller, Bionet, Mortala, Burdick		HP OEM ON:40419A	V0013C-A	Adult suction Ag/AgCl electrode, - perfect fit for <u>3.0mm Din plug</u> or <u>0.156-in.(4mm) banana plug</u> or <u>(4mm)Snap</u> leadwires - Latex-free, - Adult size replacement bulb for all 15mm diameter electrodes. Bulb volume is 5cc. - Metal Electrodes: <u>Ag/AgCl</u> . - 1 box = 1 bulb	B2
		HP OEM ON:40422A	V0013C-P	Paediatric suction Ag/AgCl electrode, - perfect fit for <u>3.0mm Din plug</u> or <u>0.156-in.(4mm) banana plug</u> or <u>(4mm)Snap</u> leadwires - Latex-free - Adult size replacement bulb for all 15mm diameter electrodes. Bulb volume is <u>2cc</u> . - Metal Electrodes: <u>Ag/AgCl</u> .	B3

● Limb Clamp ECG Electrodes,四肢夹心电电极

Compatible Brand	Product Guide	OEM P/N	MED-LINK REF NO.	Description	PRICE CODE
all Resting ECG Recommended for Nihon Konden, Fukuda, NEC, Kenz, China		Nihon kondens: H068K	V0012A	Limb clamp electrode, clothespin type - perfect fit for 3.0mm Din plug leads - Metal Electrodes: Nickel-plate. - 1 box = 4 electrodes - Color coding: IEC, red(R), yellow(L), black(RF), green(F)	A8
			V0012B-I	Limb clamp electrode, clothespin type - perfect fit for 0.156-in.(4mm) banana plug leads - Metal Electrodes: Nickel-plate.. - 1 box = 4 electrodes - Color coding:IEC, red(R), yellow(L), black(RF), green(F)	A8
			V0012D-A	Limb clamp electrode, clothespin type - perfect fit for 3.0mm Din plug leads - Metal Electrodes: Nickel-plate. - 1 box = 4 electrodes - Color coding: AHA, white(RA), black(LA), green(RL), red(LL)	A8
			V0012B	Limb clamp electrode, clothespin type - perfect fit for 0.156-in.(4mm) banana plug leadwires - Metal Electrodes: Nickel-plate. - 1 box = 4 electrodes - Color coding: AHA, white(RA), black(LA), green(RL), red(LL)	A8
all Resting ECG Recommended for GE-Marquette, Schiller, Bionet, Mortala, Burdick		HP OEM ON:40494E	V0012M-I	Adult Limb clamp Ag/AgCl electrode, clothespin type - perfect fit for 3.0mm Din plug or 0.156-in.(4mm) banana plug or (4mm)Snap leadwires - Metal Electrodes: Ag/AgCl. - 1 box = 4 electrodes - Color coding standard: IEC, red(R), yellow(L), black(RF), green(F)	B2
	same as above, but the color coding standard is IEC, white(RA), black(LA), green(RL), red(LL)		V0012M-A	Adult Limb clamp Ag/AgCl electrode, clothespin type - perfect fit for 3.0mm Din plug or 0.156-in.(4mm) banana plug or (4mm)Snap leadwires - Metal Electrodes: Ag/AgCl. - 1 box = 4 electrodes - Color coding : AHA, white(RA),black(LA),green(RL),red(LL)	B3
all Resting ECG Recommended for GE-Marquette, Schiller, Bionet, Mortala, Burdick			V0012M-I/P	Paediatric Limb clamp Ag/AgCl electrode, clothespin type - perfect fit for 3.0mm Din plug or 0.156-in.(4mm) banana plug or (4mm)Snap leadwires - Metal Electrodes: Ag/AgCl. - 1 box = 4 electrodes - Color coding standard: IEC, red(R), yellow(L), black(RF), green(F)	B2
	same as above, but the color coding standard is IEC, white(RA), black(LA), green(RL), red(LL)		V0012M-A/P	Paediatric Limb clamp Ag/AgCl electrode, clothespin type - perfect fit for 3.0mm Din plug or 0.156-in.(4mm) banana plug or (4mm)Snap leadwires - Metal Electrodes: Ag/AgCl. - 1 box = 4 electrodes - Color coding : AHA, white(RA),black(LA),green(RL),red(LL)	B3

Earth Cable

Compatible Brand	Product Guide	OEM P/N	MED-LINK REF NO.	Description	PRICE CODE
all patient monitor and cardiographs with standard Earth terminal			V0020A	Earth Cable, 6mm Jack >> Clip, - Length: 4000mm - Wires type: UL 1015 14AWG 105° 600V . - 1 bag = 1 pce - Color: yellow/green	A4
			V0020B	Earth Terminal, 6mm diameter(perfect fit for 6mm jack of V0020A) - Length: 30mm - Metal terminal: Nickel-plate. - 1 bag = 1 pce	A1