








Type Disignation

P □□□-□□□-□□



MED-LINK REF NO.	Productor Guide	Description, Name, Length, Connector type	PRICE CODE
P281-280-13	 281 280	Patient Return Plate Cable, Length: 13 ft. (4 m) PLUG: 2.35 jack	A3
P053-282-13	 053 282	Patient Return Plate Cable, Length: 13 ft. (4 m) PLUG: 6.35 Mono Plug	A3
P283-284-10	 284 283	Bipolar Electrode Cable, Length: 10 ft. (3 m) PLUG: 2*Banana Plug	A3
P285-286-09	 285 286	product name: Disposable Electrosurgical pencil, Cautery pencil; Categories: Surgical apparatus Length of handle(mm):170mm; Length of blade(mm):40mm; Diameter of blade:2.3mm; Length of wire(m): 9 ft (2.8m); Package way:100pcs/box; Size of box:46*23*28CM; Gross weight:8.0kg; Net weight:7.0kg Maximal output current:1A; Frequency confine:0-1.0MHZ; Current resistance:200; Ways of aseptic technique:Sterile eo; Term of sterility:2 year	A3
287		Disposable Electrosurgical Patient Plate-Monopolar	A2
288		Disposable Electrosurgical Patient Plate-Bipolar	A2