




Instrument Connector



Yoke type























transducers connector



















































Type Disignation

























X □ □ □ □ □

- Yoke type: A: Appott 6pin; B: Utah 4pin; C: Edward 6pin; D: BD 7pin; E: Medix 6pin F: Philips 4pin
- Compatible Brand code, as 0014 indicated Marquette, 0015 indicated Spacelabs, 0013 indicated Philips
- X: Indicated IBP Cable and Pressure Transducers












Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Plug)	INSTRUMENT CONNECTOR	Yoke (transducers side)	PRICE CODE
Chenwei			X0017A	Appott Transducer Adapter Cable, 12 ft. (3.6M), Readel 4pin>6p6c jack,	251		C0
			X0017B	Utah Transducer Adapter Cable, 12 ft. (3.6M), Readel 4pin>Din 2.0 4pin			C0
			X0017C	Edward Transducer Adapter Cable, 12 ft. (3.6M), Readel 4pin>6pin,			C0
			X0017D	BD Transducer Adapter Cable, 12 ft. (3.6M), Readel 4pin>7p jack,			E5
			X0017E	MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Readel 4pin>6pin			E5
			X0017F	Philips Transducer Adapter Cable, 12 ft. (3.6M), Readel 4pin>Round 4p			C0
Datex			X0030A	Appott Transducer Adapter Cable, 12 ft. (3.6M), Round 10J>6p6c jack,	252		C0
			X0030B	Utah Transducer Adapter Cable, 12 ft. (3.6M), Round 10J>Din 2.0 4pin			C0
			X0030C	Edward Transducer Adapter Cable, 12 ft. (3.6M), Round 10J>6pin,			C0
			X0030D	BD Transducer Adapter Cable, 12 ft. (3.6M), Round 10J>7p jack,			E5
			X0030E	MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Round 10J>6pin			E5
			X0030F	Philips Transducer Adapter Cable, 12 ft. (3.6M), Round 10J>Round 4p			C0
Fukuda			X0026A	Appott Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>6p6c jack,	253		E5
			X0026B	Utah Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>Din 2.0 4pin			E5
			X0026C	Edward Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>6pin,			E5
			X0026D	BD Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>7p jack,			E5
			X0026E	MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>6pin			E5
			X0026F	Philips Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>Round 4p			E5

Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Plug)	INSTRUMENT CONNECTOR	Yoke (transducers side)	PRICE CODE
GE-marquette	GE Solar, Dash, Eagle, Tram 		X0014A	Appott Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 11pin>6p6c jack,	254 		C0
		700078-001	X0014B	Utah Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 11p>Din2.0 4p			C0
			X0014C	Edward Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 11pin>6pin,			C0
			X0014D	BD Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 11pin>7p jack,			E5
			X0014E	MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 11pin>6pin			E5
			X0014F	Philips Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 11pin>Round 4p			C0
Lohmeier	M010 		X0019A	Appott Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 12p>6p6c jack,	079 		E5
			X0019B	Utah Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 12p>Din 2.0 4pin			E5
			X0019C	Edward Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 12p>6pin,			E5
			X0019D	BD Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 12p>7p jack,			E5
			X0019E	MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 12p>6pin			E5
			X0019F	Philips Transducer Adapter Cable, 12 ft. (3.6M), Rectangle 12p>Round 4p			E5
Mindray	PM6000,6800 		X0016A	Appott Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>6p6c jack,	255 		C0
			X0016B	Utah Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>Din 2.0 4pin			C0
			X0016C	Edward Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>6pin,			C0
			X0016D	BD Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>7p jack,			E5
			X0016E	MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>6pin			E5
			X0016F	Philips Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>Round 4p			C0








Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Plug)	INSTRUMENT CONNECTOR	Yoke (transducers side)	PRICE CODE
Nihon Konden			X0028A	Appott Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>6p6c jack,	256 		H1
			X0028B	Utah Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>Din 2.0 4pin			H1
			X0028C	Edward Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>6pin,			H1
			X0028D	BD Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>7p jack,			H1
			X0028E	MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>6pin			H1
			X0028F	Philips Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>Round 4p			H1
Philips	M1567A,M1568A kits 		X0018A	Appott Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>6p6c jack,	255 		C0
			X0018B	Utah Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>Din 2.0 4pin			C0
			X0018C	Edward Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>6pin,			C0
			X0018D	BD Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>7p jack,			E5
			X0018E	MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>6pin			E5
		M1634A	X0018F	Philips Transducer Adapter Cable, 12 ft. (3.6M), Round 12pin>Round 4p			C0
Siemens	Siemens Medical:Sirecust 700, Sirecust 630 		X0001A	Appott Transducer Adapter Cable, 12 ft. (3.6M), Round 10P>6p6c jack,	257 		C0
			X0001B	Utah Transducer Adapter Cable, 12 ft. (3.6M), Round 10P>Din 2.0 4pin			C0
			X0001C	Edward Transducer Adapter Cable, 12 ft. (3.6M), Round 10P>6pin,			C0
			X0001D	BD Transducer Adapter Cable, 12 ft. (3.6M), Round 10P>7p jack,			E5
			X0001E	MEDIX Transducer Adapter Cable, 12 ft. (3.6M), Round 10P>6pin			E5
			X0001F	Philips Transducer Adapter Cable, 12 ft. (3.6M), Round 10P>Round 4p			C0

Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Plug)	INSTRUMENT CONNECTOR	Yoke (transducers side)	PRICE CODE
Siemens Drager	Siemens/Darger:SC6002XL,7000, SC8000,SC9000XL  X0002		X0002A	Appott Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>6p6c jack,	258 		E5
			X0002B	Utah Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>Din 2.0 4pin			E5
			X0002C	Edward Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>6pin,			E5
			X0002D	BD Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>7p jack,			E5
			X0002E	MEDIX Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>6pin			E5
			X0002F	Philips Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>Round 4p			E5
Spacelabs Mindray	Reusable and non-sterile. It may be gas sterilized. (Datex-Ohmeda TC-VTK.) For use with 690-0008-00, 690-0009-00, 690-0010-00, 690-0011-00, 690-0014-01, and 690-0015-01 transducers and all Spacelabs Medical monitoring systems. Use on 90496 also requires adapter cable 700-0028-00, sold separately. Mindray: PM7000, 8000,9000,6201.  X0015		X0015A	Appott Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>6p6c jack,	259 		C0
			X0015B	Utah Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>Din 2.0 4pin			C0
			X0015C	Edward Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>6pin,			C0
		690-0021-00	X0015D	BD Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>7p jack,			E5
			X0015E	MEDIX Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>6pin			E5
			X0015F	Philips Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin>Round 4p			C0
MEK (KOREA)	 X0027		X0027A	Appott Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin 30°>6p6c jack,	004 		C0
			X0027B	Utah Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin 30°>Din 2.0 4pin			C0
			X0027C	Edward Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin 30°>6pin,			C0
		690-0021-00	X0027D	BD Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin 30°>7p jack,			E5
			X0027E	MEDIX Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin 30°>6pin			E5
			X0027F	Philips Transducer Adapter Cable, 12 ft. (3.6M), AMP 6pin 30°>Round 4p			C0

ADAPTER CABLE

Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Plug)	INSTRUMENT CONNECTOR	Yoke (transducers side)	PRICE CODE
Mennen			X0025A	Adapter cable, pressure transducer interface cable, 35 cm, Round 8p>AMP 6J	015 	270 	E5
		800-030-270	X0003A	Adapter cable, pressure transducer interface cable, 35 cm, Round 10p>2* 6J	260 	269 	F5
Siemens		13646	X0004A	Adapter cable, pressure transducer interface cable, 35 cm, Rec16p> 2*Rou 10p	261 	272 	G0
Spacelabs	For use with the 90496 and 690-0021-00.	700-0028-00	X0024A	Adapter cable, pressure transducer interface cable, 122 cm/4 ft	082 	271 	F0

DISPOSABLE PRESSURE TRANSDUCER

Compatible Brand	Productur Guide	Ttransducer Connector	MED-LINK REF NO.	Description	PRICE CODE
Appott Medex Spacelabs p/n: 690-0008-00		 262	X0020A	DTX/Plus transducer kit for patientmount application, single-use NOTE: Transducer contains built-in flush device and stopcock. Requires 690-0021-00.	C6
UTAH		 263	X0020B	DTX/Plus transducer kit for patientmount application, single-use NOTE: Transducer contains built-in flush device and stopcock. Requires 690-0021-00.	C6
Edwards Baxter		 264	X0020C	DTX/Plus transducer kit for patientmount application, single-use NOTE: Transducer contains built-in flush device and stopcock. Requires 690-0021-00.	C6
BD 贝达		 265	X0020D	DTX/Plus transducer kit for patientmount application, single-use NOTE: Transducer contains built-in flush device and stopcock. Requires 690-0021-00.	C6
MEDEX		 267	X0020E	DTX/Plus transducer kit for patientmount application, single-use NOTE: Transducer contains built-in flush device and stopcock. Requires 690-0021-00.	C6
Philips		 268	X0020F	DTX/Plus transducer kit for patientmount application, single-use NOTE: Transducer contains built-in flush device and stopcock. Requires 690-0021-00.	D5