

Instrument Connector

Electrocardiograph Connector

Plug code	MED-LINK REF NO.	Description, Name, Connector color	USE / Function	Fit Cable Diameter	PRICE CODE
034	CT0057C	Locking type DB15M, Gray, Plastic	Fukuda denshi EKG	φ6.0mm	B0
165	CT0056A	Locking type DB15M, Gray, Plastic	Fukuda ME EKG	φ6.0mm	B0
158	CT0057M	Locking type DB15M, Gray, Plastic	Motala EKG	φ6.0mm	B0
222	CT0050A	Locking type DB15M, Gray, Plastic	FMT EKG	φ6.0mm	B0
032	CT0036A	Screw type DB15M, Gray, Metal	Siemens EKG	φ6.0mm	A0
024	CT0024A	DDK Locking type 16pin, Black, Metal	Kanz EKG	φ6.0mm	B6



Patient monitor Connector for ECG

Plug code	MED-LINK REF NO.	Description, Name, Connector color	USE / Function	Fit Cable Diameter	PRICE CODE
005	CT0038A	AMP Round 6pin, Black, Plastic	CSI ECG	φ5.0mm	B2
236	CT0057G	AMP Round 5pin, Black, Plastic	Philips ECG	φ5.0mm	B2
175	CT0057H	Round 10pin, White, Plastic	Mennen ECG	φ5.0mm	B0
011	CT0057B	Round 12pin, Green, Plastic	Kontron ECG	φ5.0mm	B2
006	CT0057D	Round 8pin, Black, Plastic	NihonKonden ECG	φ5.0mm	B2
237	CT0057I	Round 11pin, Black, Plastic	NihonKonden ECG	φ5.0mm	B0
023	CT0057C	Locking type round 10J, Black, Plastic	Datex ECG	φ5.0mm	B2
018	CT0057A	Round 12pin, White, Plastic	Philips ECG	φ5.0mm	B2
170	CT0057E	Round 10pin, White, Plastic	Drager SG ECG	φ5.0mm	B0
015	CT0057F	Locking type 8pin, Red, Plastic	Bionet	φ5.0mm	B2
071	CT0014C	Locing type 7pin, Black, Plastic	Novamatrix ECG	φ5.0mm	A0
014	CT0002A	Readel 1-bit 6-pin, Green, Plastic	Criticare ECG	φ5.0mm	A5
242	CT0057J	Screw type Din-5pin, Black, Metal	Siemens EKG	φ5.0mm	B0
54	CT0057K	Stereo 6.3 Plug, black, Plastic shell	PHYSIO ECG	φ5.0mm	A2



Patient monitor Connector for SpO2

Plug code	MED-LINK REF NO.	Description, Name, Connector color	USE / Function	Fit Cable Diameter	PRICE CODE
041	CT0009A	Readel 1-bit 7-pin, Blue, Plastic	BCI SpO2	φ4.0mm	A5
211	CT0008A	Readel 2-bit 7-pin, Blue, Plastic	Goldway SpO2	φ4.0mm	A5
074	CT0001A	Readel 1-bit 5-pin, Blue, Plastic	Mindray SpO2	φ4.0mm	A5
076	CT0003A	Readel 2-bit 6-pin, Blue, Plastic	Creative SpO2	φ4.0mm	A5
044	CT0010A	FGG 1B 5-pin, Blue, metal	CSI SpO2	φ4.0mm	B4
189	CT0011A	FGG 1B 7-pin, Blue, metal	INVIVO SpO2	φ4.0mm	B4
045	CT0010C	FGG 1B 6-pin, Blue, metal	NONIN SpO2	φ4.0mm	B4
043	CT0018B	FGG 2B 10-pin, Blue, metal	MENNEN SpO2	φ4.0mm	C0
046	CT0043A	FGY 2B 8-pin, Blue, Plastic	CSI SpO2	φ4.0mm	F0
068	CT0034A	Din 8-pin, Black, Plastic	Datascope SpO2	φ4.0mm	A2
219	CT0034D	Din 7-pin, Black, Plastic	Schiller SpO2	φ4.0mm	A2
071	CT0014A	Hye 7-pin, Black, Plastic	Novamatrix SpO2	φ4.0mm	A0
072	CT0013A	hye 9-pin, Black, Plastic	Simed SpO2	φ4.0mm	A0
220	CT0015A	hye 12-pin, Black, Plastic	Spacelabs SpO2	φ4.0mm	A0
062	CT0058A	Round 12-pin, Blue, Plastic	Bruker SpO2	φ4.0mm	B2
191	CT0058B	Round 10J, Gray, Plastic	Datex SpO2	φ4.0mm	B2
128	CT0058E	Round 12-pin, Blue, Plastic	Philips SpO2	φ4.0mm	B2
047	CT0058F	Round 10-pin, Blue, Plastic	NihonKohden SpO2	φ4.0mm	B2
218	CT0058C	Round 10J, Purple, Plastic	S&W SpO2	φ4.0mm	B2
289	CT0058M	Round 6pin, Black, Plastic	Siemens SpO2	φ4.0mm	B2
003	CT0058N	Round 6pin 30°, Black, Plastic	Bionet SpO2	φ4.0mm	B2
061	CT0031A	MD6M, Blue, Plastic	Biosys SpO2	φ4.0mm	A2
217	CT0058D	RJ11 6P6C, Gray, Plastic	Palco SpO2	φ4.0mm	A2
213	CT0058H	Round 8-pin, Blue, Plastic	Mennen SpO2	φ4.0mm	B2
234	CT0035A	Karon 5J, Black, Metal	Pace tech SpO2	φ4.0mm	B0
042	CT0020A	HD14-pin, Gray, Plastic	Masimo SpO2	φ4.0mm	A5
290	CT0020C	HD14-pin, 1 bit Gray, Plastic	Nellcor SpO2	φ4.0mm	B5
116	CT0020D	HD14-pin, 1 bit, Blue, Plastic	Nellcor SpO2	φ4.0mm	B5
123	CT0058I	Locking 8P, Blue, Plastic	Ohmeda SpO2	φ4.0mm	B2
230	CT0058J	Locking 8J, Blue, Plastic	Ohmeda SpO2	φ4.0mm	B2
115	CT0058K	Quick Lock 6-pin,White, Plastic	Masimo SpO2	φ4.0mm	B2
109	CT0058C	Quick Lock 6J,White, Plastic	Masimo SpO2	φ4.0mm	B2



Patient monitor Socket

OEM P/N	MED-LINK REF NO.	Description, Name, Connector color	USE / Function	PRICE CODE
	E0001B	AMP Round 6-pin socket, Black, Plastic, Flange: 30X30cm	ECG	A3
	E0001G	AMP Round 6-pin socket, Black, Plastic, Flange: 23.5X23.5cm	ECG	A3
	E0001D	AMP Round 6-pin socket, Red, Plastic, Flange: 30X30cm	IBP	A3
	E0001F	AMP Round 6-pin socket, Red, Plastic, Flange: 23.5X23.5cm	IBP	A3
	E0001O	AMP Round 5-pin socket, Gray, Plastic, Flange: 30X30cm	ECG	A3
	E0001P	AMP Round 6-pin socket, Blue, Plastic, Flange: 30X30cm	SpO2	A3
	E0001E	Locking type 12Pin Socket, Green, Flange 28X28cm, Plastic	Kontron / ECG	A6
	CT0009B	Readel 1-bit 7J Socket, Blue, Plastic	SpO2	A5
	CT0008B	Readel 2-bit 7J Socket, Blue, Plastic	SpO2	A5
	CT0001B	Readel 1-bit 5J Socket, Blue, Plastic	SpO2	A5
	CT0003B	Readel 2-bit 6J Socket, Blue, Plastic	SpO2	A5
	CT0010D	FGG 1B 5J Socket, Blue, metal	SpO2	B3
	CT0011C	FGG 1B 7J Socket, Blue, metal	SpO2	B3
	CT0010E	FGG 1B 6J Socket, Blue, metal	SpO2	B3
	CT0018D	FGG 2B 10J Socket, Blue, metal	SpO2	C5
	CT0068A	6.3 mm Stereo socket	TEMP	A1
	CT0068B	Din1.5 pin single socket, Black	EEG	A1

E0001B



E0001G



E0001D



E0001F



E0001O



E0001P



E0001E



CT0009B



CT0010D



CT0068A



CT0068B

