


Instrument Connector




Yoke type



Disposable temp. probe connector



Probe type and Code



Specifications

Temperature Range: 0 to 50°C


Accuracy: ±0.2°C from 0 to 50°C, ±0.1°C from 25 to 45°C

© Today's quality patient care requires body temperature monitor to be an integral part of any physiological measurements in the OR, ICU, CCU, or ER. MED-LINK offers a full line of reusable, sterilizable, temperature probes to measure body temperature on your critical patients.

© Available in Rectal/Esophageal and Skin Surface Styles

© Latex-free

Type Disignation






























































Probe/Yoke type: **A**: Adult Reusable Skin-urface Probes; **B**: Adult Reusable Esophageal/Rectal Probe; **C**: Paediatric Reusable Skin-urface Probes; **D**: Paediatric Reusable Esophageal/Rectal Probe; **E**: Disposable Skin-urface Probes; **F**: Disposable Esophageal/Rectal Probes; **G**: Rectangle 2p yoke; **H**: 3.5 stereo Jack yoke, **I**: Airway Temperature Probe, **J**: DC5.5X2.5 Jack.

Compatible Brand code, as 0001 indicated YSI 400 series, 0003 Philips temperature monitors; 0004: 10K

W: Indicated Medical Temperature Probes

Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Probes,Plug)	Instrument Connector or Probe connector	Yoke (probe side) or Probe type	PRICE CODE
YSI 400 Series • Reusable temperature probes with Single Thermistors • Available in Rectal/Esophageal	Philips A3 series monitor Instrument: M3926A, M3927A, M3928A, M3929A; Marquette 7000/7010,GE Eagle, Solar, Tram, Tramscope, MAC-Lab Systems	409B	W0001A	Adult Reusable Skin-surface Probe, 10ft. (3.0 m), 25° C@2.25KΩ, φ6.3 mono plug>>φ12mm disk	053		A9
		401	W0001B	Adult Reusable Esophageal/Rectal Probe, 10ft. (3.0 m), 25°C@2.25K Ω, φ6.3 mono plug>>φ 4.7mm tip	053		A9
		427	W0001C	paediatric Reusable Skin-surface Probe, 10ft. (3.0 m), 25°C@2.25KΩ, φ6.3 mono plug>>φ5.0mm disk	053		A0
		402	W0001D	paediatric Reusable Esophageal/Rectal Probe, 10ft. (3.0 m), 25°C@2.25K Ω, φ6.3 mono plug>>φ 3.0mm tip	053		A0
		4499	W0001E	Disposable Skin-surface Probe, 0.8 m, 25°C@2.25KΩ, Rectangle 2pin>>Adhesive Pad	117		A5
		4491	W0001F	Adult Esophageal/Rectal Probe, 0.8 m, 25°C@2.25KΩ, Rectangle 2pin>>φ4.0mm tip	117		A5
		E9004ZS	W0001G	Temp. probe Interconnect Cable, 8 ft. (2.4M) , φ6.3 mono plug >> Rectangle 2pin, for use with Disposable Probes	053		A6
			W0001H	Temp. probe Interconnect Cable, 8 ft. (2.4M) , φ6.3 mono plug >> φ3.5 Stereo jack, for use with Disposable Probes	053		A6
	Allows continuous temperature monitoring of gases delivered through anesthesia and respiratory circuits. Manifold outlet included with each probe.	441A	W0001I	Airway Temperature Probe, 10 ft. (3.0M) ,25°C@2.25KΩ, φ2.5 mono plug >> Type L probe, 7.6 mm diameter tapered manifold fitting			B0
	YSI 700 Series Reusable temperature probes with Dual Thermistors	Spacelabs: 90496 Module , compatible with Spacelabs P/N: 700-0031-00 Marquette 7000/7010,GE Eagle, Solar, Tram, Tramscope, MAC-Lab Systems, compatible with GEMS P/N:9445-001	709B	W0002A	Adult Reusable Skin-surface Probe, 10ft. (3.0 m), 25° C@6K/30KΩ, φ6.3 stereo plug>>φ12mm disk	054	
		701	W0002B	Adult Reusable Esophageal/Rectal Probe, 10ft. (3.0 m), 25°C@6K/30K Ω, φ6.3 stereo plug>>φ 4.7mm tip	054		B0
			W0002C	paediatric Reusable Skin-surface Probe, 10ft. (3.0 m), 25°C@6K/30KΩ, φ6.3 stereo plug>>φ5.0mm disk	054		C5
			W0002D	paediatric Reusable Esophageal/Rectal Probe, 10ft. (3.0 m), 25°C@6K/30K Ω, φ6.3 stereo plug>>φ 3.0mm tip	054		C5

Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Probes,Plug)	Instrument Connector or Probe connector	Yoke (probe side) or Probe type	PRICE CODE
10K Series • Reusable temperature probes with Single Thermistors • Available in Rectal/Esophageal	Suitable for China OEM patient monitor: Tianrong : TR900D/E; Anke : ASC553A3, ASC553; Comon : star 8000A/B/C, star 5000,5000B/C  W0004A  W0004B  W0004C  W0004D  W0004E  W0004F		W0004A	Adult Reusable Skin-surface Probe, 10ft, (3.0 m), 25°C@10KΩ, φ6.3 mono plug>>φ12mm disk	053 		A8
			W0004B	Adult Reusable Esophageal/Rectal Probe, 10ft, (3.0 m), 25°C@10KΩ, φ6.3 mono plug>>φ4.7mm tip	053 		A8
			W0004C	paediatric Reusable Skin-surface Probe, 10ft, (3.0 m), 25°C@10KΩ, φ6.3 mono plug>>φ5.0mm disk	053 		A0
			W0004D	paediatric Reusable Esophageal/Rectal Probe, 10ft, (3.0 m), 25°C@10KΩ, φ6.3 mono plug>>φ3.0mm tip	053 		A0
			W0004E	Disposable Skin-surface Probe, 0.8 m, 25°C@10KΩ, Rectangle 2pin>>Adhesive Pad	117 		A5
			W0004F	Adult Esophageal/Rectal Probe, 0.8 m, 25°C@10KΩ, Rectangle 2pin>>φ4.0mm tip	117 		A5
		E9004ZS	W0001G	Temp. probe Interconnect Cable, 8 ft. (2.4M) , φ6.3 mono plug >> Rectangle 2pin, for use with Disposable Probes	053 		A6
			W0001H	Temp. probe Interconnect Cable, 8 ft. (2.4M) , φ6.3 mono plug >> φ3.5 Stereo jack, for use with Disposable Probes	053 		A6
		Creative Shenzhen	Creative patient monitor: PC-1000,5000,9000,UP-6000,  W0005A		W0005A	Adult Reusable Skin-surface Probe, 10ft, (3.0 m), 25°C@10KΩ, Readel 2p>>φ12mm disk	273 
	W0005B			Adult Reusable Esophageal/Rectal Probe, 10ft, (3.0 m), 25°C@10KΩ, Readel 2p>>φ4.7mm tip	273 		B5
18B20 series	Probes: DALLAS 18B20  W0006A  W0006B		W0006A	Adult Reusable Skin-surface Probe, 10ft, (3.0 m), 18B20, φ6.3 stereo plug>>φ12mm disk	054 		B0
			W0006B	Adult Reusable Esophageal/Rectal Probe, 10ft, (3.0 m), 18B20, φ6.3 stereo plug>>φ6.0mm tip	054 		B0

Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Probes,Plug)	Instrument Connector or Probe connector	Yoke (probe side) or Probe type	PRICE CODE
Philips <ul style="list-style-type: none"> • Reusable temperature probes with Single Thermistors • Available in Rectal/Esophageal 	All Philips temperature monitors that accept two-prong plugs 	21078A	W0003A	Adult Reusable Skin-urface Probe, 10ft, (3.0 m), 25°C@2.25KΩ, 2pin plug>>φ12mm disk	051 		B5
		21075A	W0003B	Adult Reusable Esophageal/Rectal Probe, 10ft, (3.0 m), 25°C@2.25KΩ, 2 pin plug>>φ4.7mm tip	051 		B5
			W0003C	paediatric Reusable Skin-urface Probe, 10ft, (3.0 m), 25°C@2.25KΩ, 2pin plug>>φ5.0mm disk	051 		B7
		21076A	W0003D	paediatric Reusable Esophageal/Rectal Probe, 10ft, (3.0 m), 25°C@2.25KΩ, 2pin plug>>φ3.0mm tip	051 		B7
		21091A	W0003E	Disposable Skin-urface Probe, 0.8 m, 25°C@2.25KΩ, 3.5 mono plug>>Adhesive Pad	112 		A5
		21090A	W0003F	Adult Esophageal/Rectal Probe, 0.8 m, 25°C@2.25KΩ, 3.5 mono plug>>φ4.0mm tip	112 		A5
			W0003G	Temp. probe Interconnect Cable, 8 ft. (2.4M) , 2 pin plug >> Rectangle 2pin, for use with Disposable Probes	051 		B3
			21082A	W0003H	Temp. probe Interconnect Cable, 8 ft. (2.4M) , 2pin plug >> φ3.5 Stereo jack, for use with Disposable Probes	051 	
S&W Artema <ul style="list-style-type: none"> • Reusable temperature probes with Single Thermistors • Available in Rectal/Esophageal 	DIASCOPE G2		W0007A	Adult Reusable Skin-urface Probe, 10ft, (3.0 m), 25°C@2.25KΩ, Round 10J>>φ12mm disk	276 		B5
			W0007B	Adult Reusable Esophageal/Rectal Probe, 10ft, (3.0 m), 25°C@2.25KΩ, Round 10J>>φ4.7mm tip	276 		B5
			W0007C	paediatric Reusable Skin-urface Probe, 10ft, (3.0 m), 25°C@2.25KΩ, Round 10J>>φ5.0mm disk	276 		B7
			W0007D	paediatric Reusable Esophageal/Rectal Probe, 10ft, (3.0 m), 25°C@2.25KΩ, Round 10J>>φ3.0mm tip	276 		B7
			W0007G	Temp. probe Interconnect Cable, 8 ft. (2.4M) , Round 10J>> Rectangle 2pin, for use with Disposable Probes	276 		B5

Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Probes,Plug)	Instrument Connector or Probe connector	Yoke (probe side) or Probe type	PRICE CODE
Siemens • Reusable temperature probes with Single Thermistors • Available in Rectal/Esophageal	Siemens medical SC 6002XL, SC7000, SC8000, SC9000XL		W0008A	Adult Reusable Skin-urface Probe, 10ft, (3.0 m), 25°C@2.25KΩ, Round 7p>>φ12mm disk			B5
			W0008B	Adult Reusable Esophageal/Rectal Probe, 10ft, (3.0m),25°C@2.25KΩ, Round 7p>>φ4.7mm tip			B5
			W0008C	paediatric Reusable Skin-urface Probe, 10ft, (3.0 m), 25°C@2.25KΩ, Round 7p>>φ5.0mm disk			B7
			W0008D	paediatric Reusable Esophageal/Rectal Probe, 10ft, (3.0m), 25°C@2.25KΩ, Round 7p>>φ3.0mm tip			B7
			W0008G	Temp. probe Interconnect Cable, 8 ft. (2.4M) , Round 7p >> Rectangle 2pin, for use with Disposable Probes			B5
25°C@2.016KΩ O.S.7929	<u>O.S.7929</u>	Ref. :09 2. 059. 6 00	W0009A	Adult Reusable Skin-urface Probe, 10ft, (3.0 m), 25°C@2.016KΩ, Screw Din3p>>φ12mm disk			B0
			W0009B	Adult Reusable Esophageal/Rectal Probe, 10ft, (3.0m),25°C@2.016KΩ, Screw Din3p>>φ4.7mm tip			B0
Temp Sensor Cut Air Shields 25°C@6.867KΩ C-100	C-100 		W0010A	Adult Reusable Skin-urface Probe, 10ft, (3.0 m), 25°C@6.867KΩ, Din3p>>φ12mm disk			B0
			W0010B	Adult Reusable Esophageal/Rectal Probe, 10ft, (3.0m),25°C@6.867KΩ, Din3p>>φ4.7mm tip			B0
15K Series • Reusable temperature probes with Single Thermistors • Available in Rectal/Esophageal	 		W0011A	Adult Reusable Skin-urface Probe, 10ft, (3.0 m), 25°C@15KΩ, Din3p>>φ12mm disk			A9
			W0011B	Adult Reusable Esophageal/Rectal Probe, 10ft, (3.0m),25°C@15KΩ, Din3p>>φ4.7mm tip			A9
50K Series • Reusable temperature probes with Single Thermistors • Available in Rectal/Esophageal	 		W0012A	Adult Reusable Skin-urface Probe, 10ft, (3.0 m), 25°C@50KΩ, Din3p>>φ12mm disk			A9
			W0012B	Adult Reusable Esophageal/Rectal Probe, 10ft, (3.0m),25°C@50KΩ, Din3p>>φ4.7mm tip			A9

