

Instrument Socket



Instrument Connector



Readel air connector

Bayonet type

Locking type

Round Connector

Rectangular Connector,

Gray Shroud Connector

HP 6.0mm diameter air plug

Yoke type



Cuff connector



Male Plastic Luer Lock

Male Plastic Screw type

Female Locking type

Male Luer Slip type

Quick disconnect Locking type

Bayonet type

Type Disignation

Y □ □ □ □ □ □

Yoke code & cuff type code: **T1**:Thigh,Single-tube; **T2**:Thigh,two-tube; **L1**: Large Adult,Single-tube; **L2**: Large Adult,two-tube; **A1**: adult,Single-tube; **A2**: adult,two-tube; **C1**: Child,Single-tube; **C2**: Child,two-tube; **I1**: Infant,Single-tube; **I2**: Infant,two-tube; **N1**: New Bo,Single-tube; **N2**: New Bo,two-tube;

Instrument Connector code. If the product is cuff, the lattice is useless, but use the "000" to express. As Y000A1

Y: Indicated Non-Invasive Blood Pressure(NIBP) Cuff & Tubes

Specifications





Blood Pressure Tubing (Air Hoses)

- Replacement tubing (air hose) for sphygmomanometers, wall manometers, aneroid gauges, and automated blood pressure monitoring systems.
- Available in Adult, Pediatric, and Neonatal styles
- 8, 12, and 24-foot lengths permit positioning of the monitor in the room as needed
- Durable tubing offers extended tubing life
- Secure connection to blood pressure cuffs to prevent air leakage

Reusable Blood Pressure Cuffs





















- Gentle TPU, strong enough for repeated inflation
- Reusable for multiple patient use
- Convenient and easy to clean
- Easy-to-use range markers and index line for proper size and placement
- Extra hook and loop for added security
- Variety of connection types to fit multiple monitoring systems
- Latex-free, PVC-free

Reusable Blood Pressure Cuffs

Compatible Brand	Product Guide	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Plug)	PRICE CODE
Philips 78352A/C, 78354A/C,M1008A/B M1176A,M1350B,M3000A M3001A (with M1598B, M1599B, tubing) Colin , Datascope - Passport, Acutor; Fukada Denshi ; Spacelabs : all; Older welch-Allyn : Models-with locking luer-type connector, Criticare .; Siemens - with bayonet type connector; Mindray , Goldway ,	One-Tube "Monitor"单管型 	M1576A	Y000T1	Reusable Blood Pressure Cuffs, Thigh ,Single-tube,Limb cir.=46~66cm● 可重复使用无创血压袖带,成人大腿,单管,大腿直径:=46~66cm	A9
		M1575A	Y000L1	Reusable Blood Pressure Cuffs, Large Adult ,Single-tube, Limb cir.:33~47cm ●可重复使用无创血压袖带,加大成人单管, 手臂直径.=33~47cm	A9
		M1574A	Y000A1	Reusable Blood Pressure Cuffs, adult ,Single-tube,Bladder 13.1x23.5cm, Limb cir.=25~35cm●袖带,成人单管,气囊13.1x23.5cm, 手臂直 径:=25~35cm	A6
		M1573A	Y000C1	Reusable Blood Pressure Cuffs, Child ,Single-tube,Bladder 9.8x18cm, Limb cir.=18~26cm●袖带,儿童单管,气囊9.8x18cm, 手臂直径:=18~26cm	A6
		M1572A	Y000I1	Reusable Blood Pressure Cuffs, Infant ,Single-tube,Bladder 7.2X13cm, Limb cir.=10~19cm●袖带,婴儿单管,气囊7.2X13cm, 手臂直径:=10~19cm	A5
		M1571A	Y000N1	Reusable Blood Pressure Cuffs, New Bo ,Single-tube,Bladder 4x9cm, Limb cir.=6~11cm●袖带,新生儿单管,气囊4x9cm, 手臂直径:=6~11cm	A5
GE-Marquette, Datex-Ohmeda GE(Marquette) with male quick-connect type connector, Dinama plus, XL,SX/XP,8100 series, Each GE corometrics; Critikon ; Datex- Ohmeda -with one male and one female quick connect type connector; CAS; New Welch-Allyn vital signs models- with threaded type Critikon style connector	Two-Tube "Monitor"双管型 	002796	Y000T2	Reusable Blood Pressure Cuffs, Thigh ,two-tube,Limb cir.=46~66cm●可 重复使用无创血压袖带,成人大腿,双管,大腿直径:=46~66cm	A9
		002791	Y000L2	Reusable Blood Pressure Cuffs, Large Adult ,two-tube, Limb cir.:33~47cm ●可重复使用无创血压袖带,加大成人双管, 手臂直径.=33~47cm	A9
		002774	Y000A2	Reusable Blood Pressure Cuffs, adult ,two-tube,Bladder 13.1x23.5cm, Limb cir.=25~35cm●袖带,成人双管,气囊13.1x23.5cm, 手臂直径:=25~35cm	A6
		002779	Y000C2	Reusable Blood Pressure Cuffs, Child ,two-tube,Bladder 9.8x18cm, Limb cir.=18~26cm●袖带,儿童双管,气囊9.8x18cm, 手臂直径:=18~26cm	A6
		002781	Y000I2	Reusable Blood Pressure Cuffs, Infant ,two-tube,Bladder 7.2X13cm, Limb cir.=10~19cm●袖带,婴儿双管,气囊7.2X13cm, 手臂直径:=10~19cm	A5
		002783	Y000N2	Reusable Blood Pressure Cuffs, New Bo ,two-tube,Bladder 4x9cm, Limb cir.=6~11cm●袖带,新生儿双管,气囊4x9cm, 手臂直径:=6~11cm	A5
Tourniquet 止血带			Y000TO	Tourniquet, LXW: 70X1100mm, 止血带, LXW: 70X1100mm	B2
Pressure Infusor 加压袋			Y000PI	For infection control purposes recommended for single patient use.	F2

Pressure Cuff Interconnect Tubing

Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Plug)	Instrument Connector	Yoke (cuff side)	PRICE CODE
Philips	Connects all sizes of Philips' reusable and disposable adult and pediatric NIBP cuffs to the monitor. Use with adult and pediatric cuffs only. CANNOT BE USED WITH NEONATAL CUFFS. Use with all Philips models except A1/A3	M1598B	YA25A06-05	Pressure Cuff Interconnect Tubing, Single tube, 5ft.1.5M, 6.0 Air plug > Bayonet type, perfect fit for 5082-184 cuff connector			B9
		M1599B	YA25A06-10	Pressure Cuff Interconnect Tubing, Single tube, 10ft. 3.0M, 6.0 Air plug > Bayonet type, perfect fit for 5082-184 cuff connector			B0
Philips	NIBP hose connects to 40401X, M45xx, M157x series pressure cuffs. Use with Philips Instrument: A1/A3, M3921A, M3923A, M3926A, M3927A, M3928A, M3929A	M3918A	YA06A06-10	Pressure Cuff Interconnect Tubing, Single tube, 10 ft 3.0 M, Bayonet type > Bayonet type perfect fit for 5082-184 cuff connector			B0
Simens Mindray Spacelabs MEK,NPB	Spacelabs:90496 Module	Spacelabs P/N: 714-0018-00	YA06A06-10	Pressure Cuff Interconnect Tubing, Single tube, 10 ft 3.0 M, Bayonet type > Bayonet type perfect fit for 5082-184 cuff connector			B0
Goldway			YA10A10-12	Blood Pressure Tubing (Air Hoses), single-tube, Adult-Pediatric, 12 ft., Screw Connector to Screw Connector,			B0
Creative Shenzhen			YA18A01-10	Pressure Cuff Interconnect Tubing, Single tube, 3.0M, Readel Air plug> Barbed to barbed connector			B0
Colin, Criticare, CSI Datascope, NihonKohden, Datascope			YA07A08-10	Pressure Cuff Interconnect Tubing, Single tube, 3.0M, Readel Air plug> Barbed to barbed connector			B0
GE Dinamap, Pro, MPS	DINAMAP, Pro, Compact, MPS Series	107365	YA24A10-12	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 12 ft., Gray Shroud to Screw Connector, Gray			C0
		107366	YA24A10-24	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 24 ft., Gray Shroud to Screw Connector,			D5
		107368	YA24A16-12	Blood Pressure Tubing (Air Hoses), two- tube, Neonatal, 12 ft., Gray Shroud to Female Slip Luer,			C0
		107363	YA24A12-12	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 12 ft., Gray Shroud to Submin Connector, Gray			C0

Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Plug)	Instrument Connector	Yoke (cuff side)	PRICE CODE
GE marquette	Eagle 3000, Dash 2000/3000/4000, Tram 400/450/451, Tram 800. 850/851. Corometrics 120	9461-217	YA23A10-12	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 12 ft., Gray, Rectangular Connector to Screw Connector,			C0
		414873-001	YA23A12-12	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 12 ft., Gray, Rectangular Connector to Submin Connector,			C0
		414874-001	YA23A16-08	Blood Pressure Tubing (Air Hoses), two- tube, Neonatal, 8 ft., Rectangular Connector to Female Slip Luer,			C0
GE medical Critikon, Kontron, Invivo	DINAMAP 1846, 8100, 9300 and Plus Series Monitors	8841	YA10A10/2-12	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 12 ft., Screw Connector to Screw Connector,			C0
		88845	YA10A10/2-24	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 24 ft., Screw Connector to Screw Connector,			D5
		8841	YA10A16-08	Blood Pressure Tubing (Air Hoses), two- tube, Neonatal, 8 ft., Screw Connector to Female Slip Luer,			C0
		88845	YA10A12-08	Blood Pressure Tubing (Air Hoses), two- tube, Neonatal, 8 ft., Screw Connector to Female Slip Luer,			C0
GE medical Corometrics GE Datex	Eagle 4000, Tram 100/200/250/300/350, Corometrics 120	9461-203	YA22A12-12	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 12 ft., Screw Connector to Submin Connector,			C0
		9461-214	YA22A10-12	Blood Pressure Tubing (Air Hoses), two- tube, Adult-Pediatric, 12 ft., Screw Connector to Screw Connector,			C0
		9461-205	YA22A16-08	Blood Pressure Tubing (Air Hoses), two- tube, Neonatal, 8 ft., Round Connector to Female Slip Luer			C0

Air Hose Connectors (Cuff side)

Air Hose Connectors

- Air hose connectors to convert your air hose to use Critikon blood pressure cuffs
- Variety of shapes and styles
- Economical way to standardize your blood pressure monitors throughout your facility.

Compatibility	Product Guide	OEM P/N	MED-LINK REF NO.	Description	Perfect fit for	PRICE CODE
Unfitted			HL0024B	Single Tubing, Black , Inner diameter: 5/32 in. Out diameter: 6mm, Purchasing Unit: M - Material: Silicone	all Barbed diameter 5/32 in. air connector	A2
Unfitted		GE P/N:300 959	SL0008A	Single Tubing, Gray , Inner diameter: 5/32 in. Out diameter: 6mm, Purchasing Unit: M - Material: TPU	all Barbed diameter 5/32 in. air connector	A2
Unfitted			SL0008B	Two-Tubing, Gray, Inner diameter: 5/32 in. Out diameter: 6mm*2, Purchasing Unit: M - Material: TPU	all Barbed diameter 5/32 in. air connector	A5
Air Hose Connector		GE P/N:300 670	A01	5/32 in. ID Barbed to 5/32 in. ID Barbed, Plastic	all Barbed diameter 5/32 in. air connector	A1
L & T		GE P/N:300 667 Welch-allyn P/N:508 2-176	A02	5/32 in. ID Barbed to Male Slip Luer, Metal	Female Locking Luer and Screw Connectors	A1
L & T		GE P/N:300 669 Welch-allyn P/N:508 2-169	A03	5/32 in. ID Barbed to Female Slip Luer, Metal	Female Locking Luer and Screw Connectors	A2
Philips (HP), Siemens, Datascope and Colin		GE P/N:330 059 Welch-allyn P/N:508 2-184	A04	Male Bayonet to 5/32 in. ID Barbed, Metal	Female Bayonet	A1
		Welch-allyn P/N:508 2-184	A05	Male Bayonet to 5/32 in. ID Barbed, Plastic	Female Bayonet	A1
		GE P/N:330 060	A06	5/32 in. ID Barbed to Female Bayonet, Metal	Male Bayonet	A4

Compatibility	Product Guide	OEM P/N	MED-LINK REF NO.	Description	Perfect fit for	PRICE CODE
Spacelabs, Datascope and Colin		GE P/N:300 668 Welch- allyn P/N:508 2-168	A07	(Plastic Luer Lock) Female Locking Luer to 5/32 in. ID Barbed, Plastic	Male Locking Luer	A1
		GE P/N:300 665 Welch- allyn P/N:508 2-165	A08	5/32 in. ID Barbed to Male Locking Luer (Plastic Luer Lock), Plastic	Female Locking Luer	A1
DINAMAP		GE P/N:300 664 Welch- allyn P/N:508 2-164	A09	Male Screw to 5/32 in. ID Barbed, Plastic	Female Screw	A1
		GE P/N:300 619 Welch- allyn P/N:508 2-161	A10	1/8 in. ID Barbed to Female Screw, Metal	Male Screw	A2
GE, Marquette and Protocol		GE P/N:330 064 Welch- allyn P/N:508 2-182	A11	(Quick Disconnect Locking type) Female Submin to 1/8 in. ID Barbed, Plastic	Male Submin	A1
GE, Marquette, Datex-Ohmeda		GE P/N:330 090	A12	(Quick Disconnect Locking type) Male Submin to 5/32 in. ID Barbed, Plastic	Female Submin	A1
Marquette one-tube systems		GE P/N:330 057 Welch- allyn P/N:508 2-181	A13	5/32 in. ID Barbed to Female Marquette locking type Connector, Plastic	male Marquette locking type Connector	A7
			A14	5/32 in. ID Barbed to male Marquette locking type Connector, Plastic	Female Marquette locking type Connector	A2
ADAPTER			A15	5/32 in. ID Barbed X2 to Female Slip Luer,, Plastic	Marquette Neonatal cuff connector	A5
			A24	5/32 in. ID Barbed X2 to 5/32 in. ID Barbed X2, Plastic	5/32 in. ID Tube	A5
Eagle 3000, Dash 2000/3000/4000, Tram 400/450/451, Tram 800/850/851			A16	5/32 in. ID Barbed to Female Slip Luer,, Plastic	Marquette Neonatal cuff connector	A5

Air Hose Connectors (Instrument side)

Compatibility	Productor Guide	OEM P/N	MED-LINK REF NO.	Description	Perfect fit for	PRICE CODE
Readel			A17	5/32 in. ID Barbed to Readel air socket, Plastic	Readel air plug	A3
			A18	5/32 in. ID Barbed to Readel air plug, Plastic	Readel air socket,	A3
Philips			A19	Male Bayonet socket to 5/32 in. ID Barbed, Metal, L=29mm	Female Bayonet	A2
			A20	Male Bayonet socket to 5/32 in. ID Barbed, Metal, L=38mm	Female Bayonet	A2
		GE P/N:330 060	A06	5/32 in. ID Barbed to Female Bayonet, Metal	Male Bayonet socket	A4
Marquette		GE P/N:330 057	A13	5/32 in. ID Barbed to Female Marquette locking type socket, Plastic	male Marquette locking type Connector	A7
			A14	5/32 in. ID Barbed to male Marquette locking type Connector, Plastic	Female Marquette locking type socket	A2
Eagle 4000, Tram 100/200/250/300/350, Corometrics 118			A22	5/32 in. ID Barbed*2 to Round Connector	Perfect fit for two-tubes	A7
Eagle 3000, Dash 2000/3000/4000, Tram 400/450/451, Tram 800/850/851, Corometrics 120			A23	Rectangular Connector,	Perfect fit for two-tubes	A7
DINAMAP, Pro, Compact, MPS Series			A24	Gray Shroud Connector,	Perfect fit for two-tubes	A7