

**Product Facts**

- Superior Data Transmission.
- Shielding, Anlistatig, Low noise.
- Tough high flexing cycle.
- Superior aging performance and enduring.
- Latex-free
- Jacket material can be medical grade PVC, TPR or TPU.
- Biocompatibility: All patient contact materials are latex free and have been tested according to ISO 10993-5,ISO 10993-10.



| MED-LINK REF NO. | Name,Description Color. Packing: 305 m / Roll (Unit price is according to meter.) | Function | PRICE CODE |
|------------------|---|----------|------------|
| DX0053A          | Single core Lead wire, (4*Tinsel)1C+TPU Jacket, cool gray, O.D.: 1.5mm,           | EEG, ECG | A2         |

TPU Jacket, Gray, OD:1.5mm

Conductor: 4\*Tinsel

| MED-LINK REF NO. | Name,Description Color, Packing: 305 m / Roll (Unit price is according to meter.) | Function | PRICE CODE |
|------------------|---|----------|------------|
| DX0012A          | Single core Lead wire,(7*Tinsel)1C+TPU Jacket, cool gray, O.D.: 2.0mm,            | ECG      | A3         |

TPU Jacket, Gray, OD:2.0mm

Conductor: 7\*Tinsel

| MED-LINK REF NO. | Name,Description Color, Packing: 305 m / Roll (Unit price is according to meter.)          | Function   | PRICE CODE |
|------------------|--|------------|------------|
| DX0110A          | Single core Shield lead wire, 28AWG/1C+Carbon film+Spiral+TPU Jacket,Cool Gray,O.D.:2.0mm  | HOLTER,ECG | A3         |
| DX0109A          | Single core Shield lead wire, 26AWG/1C+Carbon film+Spiral+TPU Jacket,Cool Gray,O.D.:2.6mm  | ECG,EKG    | A3         |
| DX0023A          | Single core Shield lead wire, 26AWG/1C+Carbon film+Spiral+TPU Jacket,Cool Gray,O.D.:3.25mm | ECG,EKG    | A3         |

TPU Jacket


Shield Spiral

Carbon Film


Conductor

Insulation


| MED-LINK REF NO. | Name,Description Color, Packing: 305 m / Roll (Unit price is according to meter.) | Function | PRICE CODE |
|------------------|---|----------|------------|
| DX0009A          | 5 cores ECG cable, 26AWG/5C+Carbon film+Spiral+TPU Jacket,Cool Gray,O.D.:5.0mm    | ECG      | A4         |



| MED-LINK REF NO. | Name,Description Color, Packing: 305 m / Roll (Unit price is according to meter.) | Function | PRICE CODE |
|------------------|---|----------|------------|
| DX0022A          | 10 cores EKG cable, #26/10C+Carbon film+Braid+PVC Jacket,Cool Gray, O.D.:6.0mm    | EKG      | A5         |



| MED-LINK REF NO. | Name,Description Color, Packing: 305 m / Roll (Unit price is according to meter.)   | Function | PRICE CODE |
|------------------|---|----------|------------|
| DX0033A          | 3+2 cores sensor cable, (#30/2C+Spiral)1C+#30/3C+Spiral+TPU Jacket,Gray, O.D.:4.0mm | SpO2     | A4         |



| MED-LINK REF NO. | Name,Description Color, Packing: 305 m / Roll (Unit price is according to meter.)         | Function | PRICE CODE |
|------------------|---|----------|------------|
| DX0047A          | 2+3Core flat sensor cable, (#30/2P+Braid)1C+#30/3C+Braid+TPU Jacket, Blue, O.D.: 2.5mmX2F | SpO2     | A4         |

