



### Connector type


















### Sensor type

S0010A Series Adult finger clip SpO <sub>2</sub> sensor Application: Adult>40kg Application site: Index or other finger	
S0022 Series Pediatric finger clip SpO <sub>2</sub> sensor Application: Pediatric>10-40kg Application site: Index or other finger	
S0012A Series Adult silicone SpO <sub>2</sub> sensor Application: Adult>40kg Application site: Index or other finger	
S0012A Series Infant silicone SpO <sub>2</sub> sensor Application: Pediatric>10-40kg Application site: Index or other finger	
S0011A Neonate silicone Wrap SpO <sub>2</sub> Sensor Application: Neonate 1-4kg, Pediatric 10-40kg; Infant 3-15kg Application site: Adult, Pediatric-Index or other finger; Neonate-foot; Infant- great toe	
S0011A Series Adult ear Clip SpO <sub>2</sub> Sensor Application: Adult>40kg Application site: ear	
S0011A Series Forehead SpO <sub>2</sub> Sensor Application: Adult>40kg Application site: Forehead	
S0011A Series Animal SpO <sub>2</sub> Sensor Application: Dog, Cat Application site: ear, tongue	

### Type Disignation



\* Masimo sensor & Nellcor oximax sensor refurbished or repair only.  
\* The brands or registered trade mark listed in the files are belong to the respected owners of the properties.

Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	OEM P/N	MED-LINK REF NO.	SENSOR TYPE	EMITTER	PRICE CODE
<b>BCI</b> (Biochem International)	<b>BCI(Biochem International):</b> 3401,3304,3303,3302,3301,3300;  <b>Spacelabs Burdick:</b> PRO2;  <b>Dräger:</b> Dräger Oxipac 2000		3044	S0076B-S	 Adult finger clip	660/905-2	C0
			3178	S0076C-S	 Pediatric finger clip	660/905-2	C0
			3044	S0076D-S	 Adult silicone soft	660/905-2	D3
			3178	S0076E-S	 Pediatric silicone	660/905-2	D5
				S0076F-S	 Neonate silicone wrap	660/905-2	D3
				S0076G-S	 Adult ear clip	660/905-2	D0
				S0076H-S	 Animal SpO2 sensor	660/905-2	D0
				S0076I-S	 Forehead SpO2 sensor	660/905-2	D0
<b>BIONET KOREA</b>	BM3			S0084B-S	Adult finger clip	660/905-2	D5
				S0084C-S	Pediatric finger clip	660/905-2	D0
				S0084D-S	Adult silicone soft	660/905-2	F5
				S0084E-S	Pediatric silicone	660/905-2	F5
				S0084F-S	Neonate silicone wrap	660/905-2	F5
				S0084G-S	Adult ear clip	660/905-2	F5
				S0084H-S	Animal SpO2 sensor	660/905-2	F5
				S0084I-S	Forehead SpO2 sensor	660/905-2	F5
<b>CSI</b> Criticare Systems International	Alaris (IVAC):VitalCheck 4400 CSI(Criticare Systems International) 503DX,503Spot,504DX,506,506DX 507ELC2 ScholarIII,8100Poet Plus, Poet (602-3A) (Sensor + Cable), Scholar (Sensor + Cable)		511SD	S0044B-S	Adult finger clip	660/940-3	D2
				S0044C-S	Pediatric finger clip	660/940-3	D2
				S0044D-S	Adult silicone soft	660/940-3	D0
				S0044E-S	Pediatric silicone	660/940-3	E5
				S0044F-S	Neonate silicone wrap	660/940-3	E0
				S0044G-S	Adult ear clip	660/940-3	F5
				S0044H-S	Animal SpO2 sensor	660/940-3	F5
				S0044I-S	Forehead SpO2 sensor	660/940-3	F5
				S0052B-S	Adult finger clip	660/940-3	D2
				S0052C-S	Pediatric finger clip	660/940-3	D2
<b>DATASCOPE</b>	Datascope:Accutorr 3/4 SAT Accutorr Plus, Passport (Database SpO2) Compatible for 0020-00-071-01			S0052D-S	Adult silicone soft	660/940-3	D0
				S0052E-S	Pediatric silicone	660/940-3	E5
				S0052F-S	Neonate silicone wrap	660/940-3	E0
				S0052G-S	Adult ear clip	660/940-3	F5
				S0052H-S	Animal SpO2 sensor	660/940-3	F5
				S0052I-S	Forehead SpO2 sensor	660/940-3	F5
				S0036B-S	Adult finger clip	660/905-2	D2
				S0036C-S	Pediatric finger clip	660/905-2	D2
<b>DATEX</b>	Datex:AS/3,CS/3,CH/S,CH/2S,CH-RS, CG/2S,Cardiocap I,Cardiocap II, Cardiocap 5, Capnomac Ultima, Light Monitor, M-Estpr, Oscar II, OSP-200, Satlite Trans, Satlite Trans Plus, System 5, ULTS00, ULTSIO, ULTSVI, ULTSVO original factory code :SAF-F4		SAF-F	S0036D-S	Adult silicone soft	660/905-2	D8
				S0036E-S	Pediatric silicone	660/905-2	D0
				S0036F-S	Neonate silicone wrap	660/905-2	D5
				S0036G-S	Adult ear clip	660/905-2	D0
				S0036H-S	Animal SpO2 sensor	660/905-2	D0
				S0036I-S	Forehead SpO2 sensor	660/905-2	D0
				S0014B-S	Adult finger clip	660/940-3	D2
				S0014C-S	Pediatric finger clip	660/940-3	D2
<b>OHMEDA</b>	Datex-Ohmeda:3740, 3710, 4700 OxiCap, 5250RGM, Biox 3760, Ohmeda 3700, 3770,3775, Rascal II), <b>GE Medical Systems:</b> 7000 Series, Corometrics 511, TRAM A.R. , TRAM Modular x00SL Series), as 200SL,400SL,600SL,800A,800SL Mennen:Horizon XL800 (Ohmeda SpO2)			S0014D-S	Adult silicone soft	660/940-3	D8
				S0014E-S	Pediatric silicone	660/940-3	D0
				S0014F-S	Neonate silicone wrap	660/940-3	D5
				S0014G-S	Adult ear clip	660/940-3	D0
				S0014H-S	Animal SpO2 sensor	660/940-3	D0
				S0014I-S	Forehead SpO2 sensor	660/940-3	D0
				S0123B-S	Adult finger clip	660/940-3	D2
				S0123C-S	Pediatric finger clip	660/940-3	F5
<b>OHMEDA</b>	Datex-Ohmeda:3740, 3710, only use for Ohmeda oximeter EMITTER: 660/940-3			S0123D-S	Adult silicone soft	660/940-3	F0
				S0123E-S	Pediatric silicone	660/940-3	G5
			OXY-W-UN	S0123F-S	Neonate silicone wrap	660/940-3	G5
			OXY-E-UN	S0123G-S	Adult ear clip	660/940-3	G5
				S0123H-S	Animal SpO2 sensor	660/940-3	G5
				S0123I-S	Forehead SpO2 sensor	660/940-3	G5

Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	OEM P/N	MED-LINK REF NO.	SENSOR TYPE	EMITTER	PRICE CODE			
<b>DOLPHIN (OSI)</b>	Dolphin ONE Oximeters	<b>131</b> 	210	S0099B-S	 Adult finger clip	660/905-2	D5			
	2150, 2100			S0099C-S	 Pediatric finger clip	660/905-2	D0			
	compatible for :2G00751.			S0099D-S	 Adult silicone soft	660/905-2	F5			
	original factory code : 201			S0099E-S	 Pediatric silicone	660/905-2	F5			
	original factory code : 360		360	S0099F-S	 Neonate silicone wrap	660/905-2	F5			
	*digital PD output directly DS.			S0099G-S	 Adult ear clip	660/905-2	F5			
				S0099H-S	 Animal SpO2 sensor	660/905-2	F5			
				S0099I-S	 Forehead SpO2 sensor	660/905-2	F5			
<b>MEK</b>	MP-100, MP110, MP400. MP500,MP600,MP1000	<b>078</b> 		S0018B-S	Adult finger clip	660/940-3	C0			
				S0018C-S	Pediatric finger clip	660/940-3	C0			
				S0018D-S	Adult silicone soft	660/940-3	D0			
				S0018E-S	Pediatric silicone	660/940-3	D0			
				S0018F-S	Neonate silicone wrap	660/940-3	D0			
				S0018G-S	Adult ear clip	660/940-3	E5			
				S0018H-S	Animal SpO2 sensor	660/940-3	E5			
<b>NELLCOR</b>	<p><b>Artema</b>: Diascope <b>Datascope</b>: Accutorr Plus,Expert,Passport;<b>Dräger</b>: Cicero,Narkomed,Vitalert 1000 <b>GEMS: Corometrics</b> 118, 120, 129, 188, 510, 556; <b>GEMS: Critikon</b> 9710, 9720 Plus, Compact, Dinamap Select ,MPS, Oxyshuttle, Pro 300/400 (non-OxiMax®); <b>GEMS</b>: Eagle Series, Dash Series, Hellige Eagle 1000, 3000, Eagle SMK ,SMU, SMV, SOLAR Series, TRAM Modular x50SL Series; <b>Invivo Research</b>: M12 (Nellcor® SpO2), Millennia 4500 Plus 2; <b>Invivo/M.D.E.</b>: Escort 100 OPT11A, Escort 300A, Escort II OPT30, Prism (Nellcor® SpO2); <b>Ivy Biomedical</b>: 405A, 405D, 405P, 405T, Vital-Guard 450 CN ; <b>Mennen</b>: Horizon XL (Nellcor® SpO2) ,MR 1300 Mercury 1200; Nellcor®: N100 (req. preamp cable), N180, N200 (req. preamp cable), N-20, NPB-40, NPB-75, NPB-290/5 ,N-390, N-395, N-3000, N-6000; <b>Nihon Kohden</b>: BSM-1102 LifeScope EC, BSM-2304 iPro (LifeScope i), BSM-4102 Procyon LT (LifeScope P), BSM-4104 Procyon (LifeScope P), BSM-4112 Procyon LT (LifeScope P), BSM-4114 Procyon (LifeScope P), BSM-9510 LifeScope M (Nellcor® SpO2), OGS-2000 PocketCare (Nellcor® SpO2); Pace Tech: Vitalmax 800 Series (black face, Nellcor® SpO2);</p> <p><b>Philips</b>: 78352C, 78354C, 78834C, M1020A, M1025B, M1205A, M1350B, M1722A/B, M1732A/B, M2475B, M2601A, M3000A, M3500B, M4735A ,CMS ,Merlin Systems, V24, V26, Viridia A1, Viridia A3; <b>Siemens</b>: SC6000, SC7000, SC8000, SC9000XL, Sirecust 630 , Sirecust 700, Sirem SpO2 Module; <b>SpaceLabs</b>: 90351-0/6, 90465 (Nellcor® SpO2), 90466 (Nellcor® SpO2) ,90467 (Nellcor® SpO2), 90489 (Nellcor® SpO2) ,90496/90369/90367 Ultraview (Nellcor® SpO2) ,90651A-08 ,IM77; <b>Welch Allyn Protocol</b>: Atlas (Nellcor® SpO2), <b>Propaq CS</b> (Nellcor® SpO2); Propaq 106-EL, Propaq 202-EL, Propaq 204-EL, Propaq 206-EL, Micropaq (Nellcor® SpO2) , Spot Vital Signs, Vital Signs (non-OxiMax®)</p>	<b>048</b>  	DS-100A	S0010B-S	Adult finger clip	660/905-2	C0			
			DS-YSYSPD	S0010C-S	Pediatric finger clip	660/905-2	C0			
				S0010D-S	Adult silicone soft	660/905-2	D5			
				S0010E-S	Pediatric silicone	660/905-2	D5			
			D-YS	S0010F-S	Neonate silicone wrap	660/905-2	D5			
			D-YSE	S0010G-S	Adult ear clip	660/905-2	E5			
				S0010H-S	Animal SpO2 sensor	660/905-2	E5			
				S0010I-S	Forehead SpO2 sensor	660/905-2	E5			
			<b>NELLCOR (OXI TECH)</b>	All Nellcor & Nellcor equipped OEM monitors not listed above	<b>049</b> 	DS-100A	S0026B-S	Adult finger clip	660/905-2	E2
							S0026C-S	Pediatric finger clip	660/905-2	E3
S0026D-S	Adult silicone soft	660/905-2					E6			
S0026E-S	Pediatric silicone	660/905-2					E0			
S0026F-S	Neonate silicone wrap	660/905-2					E5			
S0026G-S	Adult ear clip	660/905-2					F0			
S0026H-S	Animal SpO2 sensor	660/905-2					F0			
	S0026I-S	Forehead SpO2 sensor	660/905-2	F0						

\* refurbished or repair only.

Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	OEM P/N	MED-LINK REF NO.	SENSOR TYPE	EMITTER	PRICE CODE
<b>NIHON KOHDEN</b>	Bedside Monitor: BSM-4101,4103  4111,4113 Life scope P  BSM-5105,5135 Life scope A  BSM-2351,2353 Life scope L  BSM-2301,2303 Life scope I	<b>067</b> 	TL-101T	S0005B-S	 Adult finger clip	660/940-3	D0
				S0005C-S	 Pediatric finger clip	660/940-3	D0
				S0005D-S	 Adult silicone soft	660/940-3	E5
				S0005E-S	 Pediatric silicone	660/940-3	E0
				S0005F-S	 Neonate silicone wrap	660/940-3	E5
				S0005G-S	 Adult ear clip	660/940-3	F0
				S0005H-S	 Animal SpO2 sensor	660/940-3	F0
				S0005I-S	 Forehead SpO2 sensor	660/940-3	F0
<b>NONIN</b>	NONIN: 8500,8600,8700,8800	<b>065</b> 	8000AA	S0080B-S	Adult finger clip	660/905-2	C0
			8000AP	S0080C-S	Pediatric finger clip	660/905-2	C0
			8000AA	S0080D-S	Adult silicone soft	660/905-2	D5
			8000AP	S0080E-S	Pediatric silicone	660/905-2	D5
				S0080F-S	Neonate silicone wrap	660/905-2	D5
				S0080G-S	Adult ear clip	660/905-2	E5
				S0080H-S	Animal SpO2 sensor	660/905-2	E5
				S0080I-S	Forehead SpO2 sensor	660/905-2	E5
<b>NOVAMETRIX</b>	Novamatrix:505,510,511,515B/C,520 7000, 7100, Oxypleth ; original match : 8776-00  Spacelabs:90351-D/G ,90465, 90466, 90467, 90489, 90651-07/08 , 90419-01, 90651A-07 original match : 8776-00	<b>226</b> 		S0062B-S	Adult finger clip	660/940-3	D5
				S0062C-S	Pediatric finger clip	660/940-3	D0
				S0062D-S	Adult silicone soft	660/940-3	F5
				S0062E-S	Pediatric silicone	660/940-3	F5
				S0062F-S	Neonate silicone wrap	660/940-3	F5
				S0062G-S	Adult ear clip	660/940-3	F5
				S0062H-S	Animal SpO2 sensor	660/940-3	F5
				S0062I-S	Forehead SpO2 sensor	660/940-3	F5
<b>Masimo (original)</b>	* original or refurbished or repair only.	<b>115</b> 	LNOP DCI	S0020B-S	Adult finger clip	660/905-2	
			LNOP DCIP	S0020C-S	Pediatric finger clip	660/905-2	
				S0020D-S	Adult silicone soft	660/905-2	
				S0020E-S	Pediatric silicone	660/905-2	
			LNOP YI	S0020F-S	Neonate silicone wrap	660/905-2	
			LNOP TC-I	S0020G-S	Adult ear clip	660/905-2	
				S0020H-S	Animal SpO2 sensor	660/905-2	
			LNOP TF-I	S0020I-S	Forehead SpO2 sensor	660/905-2	
<b>PALCO</b>	Palco: ALL	<b>217</b> 		S0079B-S	Adult finger clip	660/905-2	F5
				S0079C-S	Pediatric finger clip	660/905-2	F0
				S0079D-S	Adult silicone soft	660/905-2	H5
				S0079E-S	Pediatric silicone	660/905-2	H5
				S0079F-S	Neonate silicone wrap	660/905-2	H5
				S0079G-S	Adult ear clip	660/905-2	H5
				S0079H-S	Animal SpO2 sensor	660/905-2	H5
	S0079I-S	Forehead SpO2 sensor	660/905-2	H5			
<b>PHILIPS</b>	M3,M4,MP20,MP30,MP40,MP50, MP70,MP80,VM4,VM6,VM8	<b>052</b> 	MI191A	S0015B-S	Adult finger clip	660/905-2	D8
			MI192A	S0015C-S	Pediatric finger clip	660/905-2	D8
			MI191A	S0015D-S	Adult silicone soft	660/905-2	E5
			MI192A	S0015E-S	Pediatric silicone	660/905-2	E0
			MI193A	S0015F-S	Neonate silicone wrap	660/905-2	D0
			MI194A	S0015G-S	Adult ear clip	660/905-2	F0
				S0015H-S	Animal SpO2 sensor	660/905-2	F0
	S0015I-S	Forehead SpO2 sensor	660/905-2	F0			