




Instrument Connector

























Connector(Sensor side)



















































Type Disignation















































- The MED-LINK product range will assist you by standardizing your EXT- Cable & Adapter cable selection while providing very durable products compatible with the many leading manufacturers of patient monitor systems.
- Cable length: 4 ft ~ 8 ft
- Storage and transport environment: -5~+40°C; Operating temperature: 0~+40°C
- Biocompatibility: All patient contact materials are latex free and have been tested according to ISO 10993-5,ISO 10993-10.
- * **Masimo SpO₂ Adapter cable refurbished or repair only.**
- * **The brands or registered trade mark listed in the files are belong to the respected owners of the properties.**























Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	CONNECTOR (SENSOR SIDE)	PRICE CODE
BCI	(Biochem International): 3401,3304,3303,3302,3301,3300;		S0076BC-L	Exi-Cable, 8ft, DB9M>DB9F, Suitable for BCI probe			B4
	(Biochem International): 6100, 9100, Advisor, Autocorr, Mini Torr BCI (Biochem International), Mini Torr Plus		S0111BC-L	Exi-Cable, 8ft, Readel 1bit. 7P>DB9F, Suitable for BCI probe			B4
Biocare (China)	BM9000		S0022BC-L	Exi-Cable, 8ft, Readel 2bit. 6P>DB9F, Suitable for BCI probe			B4
Bionet (Korea)	BM3Vet, BM31		S0084BC-L	Exi-Cable, 8ft, Readel 1bit. 7P>DB9F, Suitable for Bionet probe			B4
Biosys (Korea)	BPM300		S0065NE-L	Exi-Cable, 8ft, MD6M>DB9F, Suitable for Nellcor probe			C0
Bruker (ODAM)	Bruker (Odam) all		S0071BC-L	Exi-Cable, 8ft, Round12P>DB9F, Suitable for BCI probe			C5
Colin	Datascope:Accutorr Plus, Nellcor:NPB-290/5, N-390,N-395,N-3000 N-6000,NPB-290/5,N-390,N-395,N-3000,N-6000, Nihon Kohden:BSM-4104A, BSM-4114A, Colin BP88/306	SCP-10	S0008CO-L	Exi-Cable, 8ft, AMP14P>DB9F, Suitable for Nellcor probe			D2
Creative (China)	Creative: PC9000,UP6000, TOB-MED同科博越: TOM MS-6-P, EDAN(China)理邦: M9, BIOCARE: BM9000, Goldway, CONTEC(康泰)		S0063BC-L	Exi-Cable, 8ft, Readel 2bit. 6P>DB9F, Suitable for BCI probe			B4
CSI	CSI (Criticare Systems International):506 Spot CSI ,507SD(Criticare Systems International):503 Spot		S0044BC-L	Exi-Cable, 8ft, DB9M>DB9F, Suitable for CSI probe			B3
	CSI (Criticare Systems International):Scholar , Poet (602-3A), 1100 , 581, 508, 507, 507SD,504, 503, 501; Mindray 500A		S0100BC-L	Exi-Cable, 8ft, Lemo 1B.5>DB9F, Suitable for CSI probe			D3
Datascope	Datascope:Accusat Datascope:Accutorr 3/4 SAT Datascope:Accutorr Plus (Datascope Sp02) Datascope:Passport (Database Sp02)		S0052BC-L	Exi-Cable, 8ft, DIN8P>DB9F, Suitable for Datascope probe			C5























OEM	OEM MODELS NUMBER	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	CONNECTOR (SENSOR SIDE)	PRICE CODE
DATEX	Datex-Ohmeda:AS/3, CS/3, CH/S, CH/2S, CH-RS, CG/2S, Cardiocap I, Cardiocap II, Cardiocap 5, Capnomac Ultima, Light Monitor, M-Estpr, Oscar II, OSP-200, Satellite Trans, Satellite Trans Plus , ULTSVI . ULTSVO, ULTSIO, ULTS00	OXY-C3	S0036BC-L	Exi-Cable, 8ft, Round 10J>DB9F, Suitable for DATEX probe			C4
		OXY-SL3	S0036OH-L	Exi-Cable, 8ft, Round 10J>8J, Suitable for DATEX probe			D0
Datex-Ohmeda	Suitable for Datex-Ohmeda:Ohmeda 3700, Ohmeda 3710, Ohmeda 3770, Ohmeda 3775, Biox 3740, Biox 3760, 3800 OxiCap, 4700 OxiCap, 5250RGM, Rascal II ; GE Medical Systems:Corometrics 511; Mennen:Horizon XL (Ohmeda SpO2)		S0014BC-L	Exi-Cable, 8ft, Hyc7P>DB9F, Suitable for Ohmeda probe			C22
		OXY-OL3	S0014OH-L	Exi-Cable, 8ft, Hyc7P>8J, Suitable for Ohmeda probe			D0
GE-Ohmeda	Sensor solutions for GE Ohmeda TruSat™ pulse oximeter:	OXY-MC3	S0127OH-L	Interconnect cable, 8ft DB9M>8J, Suitable for Ohmeda probe			F0
Dolphin Medical (OSI)	Dolphin ONE Oximeters	110	S0099BC-L	Exi-Cable, 8ft, AMP14P>DB9F, Suitable for Dolphin probe			C0
EDAN (China)	EDAN: M9, WANDONG 万东: BYS-8000, MINDRAY迈瑞: PM-8000		S0040BC-L	Exi-Cable, 8ft, Readel 2bit. 6P>DB9F, Suitable for BCI probe			B4
GE Medical (Nellcor module) (NELLCOR)	GE Medical Systems:Eagle GE Medical Systems:Dash GE Medical Systems:SOLAR GE Medical Systems:TRAM Modular x50SL Series	E9004GE	S0002NE-L	Exi-Cable, 8ft, Rectangle 11P>DB9F, Suitable for Nellcor probe			B5
GE Medical (Oxi-tech module)	GE Medical Information technology: Oximax module		S0122OX-L	Exi-Cable, 8ft, Rectangle 11P>DB9F, Suitable for Nellcor oximax probe			F0
GE Medical (Masim module) (Masimo)	GE Medical Systems:Dash,SOLAR Series (Masimo module)		S0025MA-L	Exi-Cable, 8ft, Rectangle 11P>Masimo 6J, Suitable for Masimo LNOP probe			D5
			S0025CO-L	Exi-Cable, 8ft, Rectangle 11P>DB9F, Suitable for Nellcor probe			D5
GE-Ohmeda	Sensor solutions for Datex-Ohmeda S/5 FM ports:	OXY-ES3	S0002OH-L	Interconnect cable, 8ft, Rectangle 11P>8J, Suitable for Ohmeda probe			D0



















Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	CONNECTOR (SENSOR SIDE)	PRICE CODE
Goldway	UT4000A, UT4000C, UT4000E, UT4000Apro, UT4000Fpro, UT6000A		S0038BC-L	Exi-Cable, 8ft, Lemo 1B. 5>DB9F, Suitable for BCI probe	044 		C5
Goldway	Goldway: UT4000F, PACE TECH		S0007BC-L	Exi-Cable, 8ft, Readel 1bit. 5p>DB9F, Suitable for BCI probe	014 		B4
INVIVO (Nellcor module) (NELLCOR模块)	4500,4500Plus,Millenia,Omega14 45,Omni-TrakTVS,Sat Scout,original cable number 9383		S0083NE-L	Exi-Cable, 8ft, Lemo 1B. 7p>DB9F, Suitable for Nellcor probe	235 		C4
INVIVO (Masim module) (Masim模块)	M8		S0094MA-L	Exi-Cable, 8ft, Lemo 1B. 7p>Masimo 6J, Suitable for Masimo LNOP probe	235 		D0
			S0094CO-L	Exi-Cable, 8ft, Lemo 1B. 7p>DB9F, Suitable for Nellcor probe	235 		D0
Kontron	Kontron:7840 Kontron:7845 Kontron:Kolormon Kontron:Micromon Kontron:Minimon 7137 Kontron:Supermon Kontron:Trakmon		S0019NE-L	Exi-Cable, 8ft, Round 12P>DB9F, Suitable for Nellcor probe	070 		C2
Masim	Masim all MDE:EscortII OPT30:Prism ,GE Marquette:,Pro1000,V300 , INVIVO:T8 , Schiller: , ZOLL: , Datascope: Trio,Passport2, Spectrum, * Original or refurbished or repair only.	PC-08	S0020MA-L	Exi-Cable, 8ft, Amp14p>Masimo 6J, Suitable for Masimo LNOP probe (ORIGINAL)	042 		
			S0020CO-L	Exi-Cable, 8ft, Amp14P>DB9F, Suitable for Nellcor probe	042 		C5
MEK (Korea)	MP-100, MP110, MP400, MP500,MP600,MP1000		S0018BC-L	Exi-Cable, 8ft, DB9M>DB9F, Suitable for MEK probe	078 		B5
	MEK: NEW		S0115BC-L	Exi-Cable, 8ft, Round 8p>DB9F, Suitable for MEK probe	227 		D0
Mennen			S0103NE-L	Exi-Cable, 8ft, Round 8p>DB9F, Suitable for Nellcor probe	228 		D0
	Horizon2000		S0009NE-L	Exi-Cable, 8ft, Lemo 2b. 8P>DB9F, Suitable for Nellcor probe	043 		F0

Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	CONNECTOR (SENSOR SIDE)	PRICE CODE
Mindray (csi module)	MEC-509,509+,503,507, Marketing:SFDA, 销售区域:中国SFDA OEM No.512A-30-06074	512A-30-06074	S0011NE-L	Exi-Cable, 8ft, Lemo 1B.5p>DB9F, Suitable for CSI probe			C5
Mindray (mindray module)	PM-7000,6201,8000,9000, Marketing:SFDA,CE,FDA; 销售区域:国内,CE,FDA, OEM No.0010-20-42594	0010-20-42594	S0040NE-L	Exi-Cable, 8ft, Readel 2bit 6p>DB9F, Suitable for Nellcor probe			B4
	PM5000,PM6000,PM8100 Marketing:SFDA,CE. 销售区域:中国,CE OEM No.0010-21-11957	0010-21-11957	S0045OX-L	Exi-Cable, 8ft, Round 12p>DB9F, Suitable for Nellcor probe			B0
	PM-6800 Marketing:SFDA,CE.FDA 销售区域:中国,CE,FDA OEM No.0010-30-42737	0010-30-42737	S0125CO-L	Exi-Cable, 8ft, Round 7p>DB9F, Suitable for Nellcor probe			F0
Mindray (masim module)	VS800,PM80002,9201,(PM-6201,7000,8000,M1K0,M2K0 since 2006), Marketing:SFDA,CE.FDA 销售区域:中国,CE,FDA OEM No.0010-30-42625	0010-30-42625	S0116CO-L	Exi-Cable, 8ft, Readel 2bit 60° 6p>DB9F, Suitable for Masimo probe			D0
	PM-9000,9100,9200,(PM-6201,7000,8000,M1K0,M2K0 before 2006) Marketing:SFDA,CE.FDA 销售区域:中国,CE,FDA OEM No.0010-30-42602	0010-30-42602	S0013CO-L	Exi-Cable, 8ft, Readel 2bit 40° 6p>DB9F, Suitable for Masimo probe			D0
	PM-6000 (masimo module, 12pin socket) Marketing:SFDA,CE.FDA 销售区域:中国,CE,FDA OEM No.0010-30-12452	0010-30-12452	S0012CO-L	Exi-Cable, 8ft, Round 12p>DB9F, Suitable for Masimo probe			D0
	PM-6800 (masimo module, 7pin socket) Marketing:SFDA,CE.FDA 销售区域:中国,CE,FDA OEM No.0010-30-42738	0010-30-42738	S00126CO-L	Exi-Cable, 8ft, Round 7p>DB9F, Suitable for Masimo probe			F0
	PM9000express, PM9000, 9100, 9200 (masimo module) Marketing:SFDA,CE.FDA 销售区域:中国,CE,FDA OEM No.9200-30-10707	9200-30-10707	S0013MA-L	Exi-Cable, 8ft, Readel 2bit 6p>DB9F, Suitable for Masimo probe			F5
	PM6000 (masimo module) Marketing:SFDA,CE.FDA 销售区域:中国,CE,FDA OEM No.6200-30-11721	6200-30-11721	S0012MA-L	Exi-Cable, 8ft, Round 12p>DB9F, Suitable for Masimo probe			F5
	Nellcor	Datascope:Accutorr Plus, Nellcor:NPB-290/5, N-390,N-395,N-3000 N-6000,NPB-290/5,N-390,N-395,N-3000,N-6000, Nihon Kohden:BSM-4104A, BSM-4114A, Colin BP88/306	SCP-10	S0001NE-L	Exi-Cable, 8ft, AMP14P>DB9F, Suitable for Nellcor probe		

Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	CONNECTOR (SENSOR SIDE)	PRICE CODE
Nellcor	Datascope: Expert,Passport, GE Medical: Corometrics 188/129,CritikonDinamap9700/9710/9720 Plus/Compact/MPS/Pro 300/Pro 400/Select, IvyBiomedical: 405P,405T, M.D.E.: Escort II OPT30,Prism, Nellcor: N100,N200,N180,N-20,NPB-40,NPB-75 Nonin: 8500,8604, Welch Allyn Protocol: Propaq 206-EL,Propaq 204-EL,Propaq 202-EL,Propaq 106EL	EC-4	S0010NE-S	Exi-Cable, 4ft, DB9M>DB9F, Suitable for Nellcor probe	059 		A8
		EC-8	S0010NE-L	Exi-Cable, 8ft, DB9M>DB9F, Suitable for Nellcor probe	059 		A8
Nellcor (Oxi-tech)	All Nellcor & Nellcor equipped OEM monitors not listed above		S0026NE-A	Adapter-Cable, 1ft, DB9M>DB9F, Suitable for Nellcor probe	240 		C0
	All Nellcor & Nellcor equipped OEM monitors not listed above	DEC-4	S0026OX-S	Exi-Cable, 4ft, DB9M>DB9F, Suitable for Nellcor oximax probe	240 		C7
		DEC-8	S0026OX-L	Exi-Cable, 8ft, DB9M>DB9F, Suitable for Nellcor oximax probe	240 		C8
	Nellcor: N550, N-595 GE: Dinamap PRO (w/Oxi) , Schiller	DOC-10	S0117OX-L	Exi-Cable, 8ft, AMP14P>DB9F, Suitable for Nellcor probe	116 		C5
Neusoft (China)	NSC series patient monitor		S0032OX-L	Exi-Cable, 8ft, Round 12p>DB9F, Suitable for Nellcor probe	128 		C0
Nihon Konden	AL-800PA Module,BSM-1101 LifeScope EC,BSM-2101/2 LifeScope L,BSM-3101 LifeScope LC,BSM-7100/7200 LifeScope 8,BSM-8200/1 LifeScope 12, BSM-8301/2 LifeScope 9,SM-8800 LifeScope14,OLV-1100A		S0069NE-L	Exi-Cable, 8ft, Round 10p>DB9F, Suitable for Nihon Konden probe	047 		C0
			S0069NI-L	Exi-Cable, 8ft, Round 10p>DB9F, Suitable for Nihon Konden probe	047 		G0
	Bedside Monitor: BSM-4101,4103,4111,4113 Life scope P,BSM-5105,5135 Life scope A,BSM-2351,2353 Life scope L,BSM-2301,2303 Life scope I	TL-201T	S0005NI-L	Exi-Cable, 8ft, Rectangle 14p>DB9F, Suitable for Nihon Konden probe	214 		F0
Nonin	8604,8604D		S0080NE-L	Exi-Cable, 8ft, Lemo 1b 6p>DB9F, Suitable for Nonin probe	045 		C5

OEM	OEM MODELS NUMBER	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	CONNECTOR (SENSOR SIDE)	PRICE CODE
Nonin	8600 series, 8500,8700,8800,9600,9700		S0109NE-L	Exi-Cable, 8ft, AMP 8p>DB9F, Suitable for Nonin probe	 215		D0
Novamatrix	Novamatrix:505, 510,511,515B/C, 520, 7000, 7100, Oxypleth ; Spacelabs Medical:90351-D/G , 90465, , 90466, 90467, 90489, 90651-07/08 , 90419-01, 90651A- 07	8776-00	S0062NE-L	Exi-Cable, 8ft, Hypertronics 7p>DB9F, Suitable for Nonin probe	 071		C2
Novamatrix	Oxypleth (Extension Cable),7100 (Extension Cable),7000 (Extension Cable),520 (Extension Cable),515B/C (Extension Cable) ,505 (Extension Cable),500		S0062HY-L	Exi-Cable, 8ft, Hypertronics 7p>Hys 7J, Suitable for Novamatrix probe	 071		D0
Pace Tech Generra	Pace Tech:All w/ black face (except Vitalmax 800 w/Nellcor® SpO2), Generra: 710. 200. 300. 400. 810, 530	4510	S0082BC-L	Exi-Cable, 8ft, Switchcraft 5J>DB9F, Suitable for BCI probe	 234		C5
			S0091BC-L	Exi-Cable, 8ft, Readel 1bit 7p>DB9F, Suitable for BCI probe	 041		B4
Philips Agilent	Philips Medical:M2601A,M3000A,M3500 B,M4735A	M1943A	S0015OX-L	Exi-Cable, 8ft, hemicycle 8p>DB9F, Suitable for Nellcor probe	 052		B0
	Philips Medical:M2601A,M3000A,M3500 B,M4735A	M1941A	S0015PH-L	Exi-Cable, 8ft, hemicycle8P>hemicycle 8J, Suitable for Philips probe	 052		C4
	Philips Medical:78352C,78354C,78834C, M1020A,M1025B,M1205A,M1350 B,M1722A/B,M1732A/B,M2475B, Viridia A1, Viridia A3	M1900B	S0003OX-L	Exi-Cable, 8ft, Round 12p>DB9F, Suitable for Nellcor probe	 128		B0
	Philips Medical:M2475B, M1732A/B, M1722A/B, M1020A, 78833C, 78833A, 78834A, 78352C, 78352A	M1940A	S0003PH-L	Exi-Cable, 8ft, Round 12P>hemicycle 8J, Suitable for Philips probe	 128		C2
S&W (Artema)	S&W(Diascope G2) Patient monitor, Diascope 1,2,3	SL033057	S0023OX-L	Exi-Cable, 8ft, Round 10J>DB9F, Suitable for Nellcor probe	 218		C5
Schiller (nellcor module)	Schiller Argus TM-7		S0118NE-L	Exi-Cable, 8ft, Din7P>DB9F, Suitable for Nellcor probe	 217		C2

Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	CONNECTOR (SENSOR SIDE)	PRICE CODE
Schiller (masim module)	Argus Pro, Argus LCM		S0028MA-L	Exi-Cable, 8ft, DIN7P>Masimo 6J, Suitable for Masimo probe			D0
SIMED (Baxter)	Baxter:ASAT , SIMED:S-50 SiMED:S-100, SIMED:S-100e	600015	S0120NE-L	Exi-Cable, 8ft, Hypertronics 8p>DB9F, Suitable for Nellcor probe			C0
Siemens Drager medical	Siemens medical SC 6002XL, SC7000, SC8000, SC9000XL		S0072BC-S	Exi-Cable, 4ft, round 7P>DB9F, Suitable for Nellcor probe			C2
	Siemens medical SC 6002XL, SC7000, SC8000, SC9000XL		S0072BC-L	Exi-Cable, 8ft, round 7P>DB9F, Suitable for Nellcor probe			C2
	Drager medical: Vitalert 1000, Cicero		S0119NE-L	Exi-Cable, 8ft, Lemo 2b 8P>DB9F, Suitable for Nellcor probe			H0
	Siemens Medical:Sirecust 700 Siemens Medical:Sirecust 630 Siemens Medical:Sirem SpO2 Module		S0078NE-L	Exi-Cable, 8ft, hemicycle10P>DB9F, Suitable for Nellcor probe			C5
Spacelabs	Spacelabs Medical:90351-0/6; 90465; 0466; 90489; 90651A-08; 90651-A; IM77,90308, PC Scout	700-0020-0	S0035NE-L	Exi-Cable, 8ft, Hypertronics 12p>DB9F, Suitable for Nellcor probe			C2
	Spacelabs Medical: 1050,1600,1700,90369,90469,90496 Ultraview	700-0030-00	S0004NE-L	Exi-Cable, 8ft, Rectangle 10P>DB9F, Suitable for Nellcor probe			C0
	Single Adapter Cable: 90351-D/G,90465,90466,90467,90489,90651-07/08 ; 2nd Adapter Cable: 90351-D/G,90465,90466,90467,90489,90651-07/08		S0035HY-L	Exi-Cable, 8ft, Hypertronics 12p>Hypertronics 7J, Suitable for Novamatrix probe		 232 Code:	D0
	Single Adapter Cable: 90351-D/G,90465,90466,90467,90489,90651-07/08 ; 2nd Adapter Cable: 90351-D/G,90465,90466,90467,90489,90651-07/08		S0035HY-A	Exi-Cable, 1ft, Hypertronics 12p>Hypertronics 7J, Suitable for Novamatrix probe		 232 Code:	D0
	Spacelabs Medical:90496 Ultraview (Novamatrix® SpO2) (Hypertronics to 10-pin)		S0004HY-A	Exi-Cable, 1ft, Rectangle10P>Hys 7J, Suitable for Novamatrix probe		 232 Code:	F0

Compatible Brand	Compatible Models No.	OEM P/N	MED-LINK REF NO.	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	CONNECTOR (SENSOR SIDE)	PRICE CODE
AnKe	AnKe安科:553; Biolight宝莱特:Cheersails千帆:CSM-2000;Goldway金科威:UT4000; WanDong万东;Jumbor精博:JB9101/9105/9107; TianRong天荣; BaoZhong保中生物; HongRunDa宏润达:EM-6; QiLan启蓝:M-III; JERRY洁瑞:PM2000; JoyceTracker今世创科: BMD-50; Mindray迈瑞:PM9000; Richer荣创:RC9000; ShenMai深迈:KP-900; Spurt世博达:CP-1200; SHIPAI狮牌: KMI-1000; ZhongHengSheng中恒生:KD580; ZhongYuan中原: HD-2000		S0007BC-L	Exi-Cable, 8ft, Readel 1bt 5P>DB9F, Suitable for BCI probe	074 		B4
Biolight	M6,M12		S0029BC-L	Exi-Cable, 8ft, Readel 1bt 5P>DB9F, Suitable for BCI probe	074 		B4
	BLT:M69, M66		S0101BC-L	Exi-Cable, 8ft, Lemo 1b 5P>DB9F, Suitable for BCI probe	044 		C5
BIONET	BM3,BM3 Plus.		S0121BC-L	Exi-Cable, 8ft, AMP 6P>DB9F, Suitable for Bionet probe	004 		C5
Choice oximax module	NELLCOR OXIMAX module		S0112OX-L	Exi-Cable, 8ft, Readel 9bt 5P>DB9F, Suitable for BCI probe	074 		C0
FMT			S0089BC-L	Exi-Cable, 8ft, Readel 2bt 6P>DB9F, Suitable for BCI probe	076 		B0
GOLDWAY	4000B		S0061BC-L	Exi-Cable, 8ft, Readel 2bt 7P>DB9F, Suitable for BCI probe	211 		B5
Hellige	Servomed 4.Generation SCP 844,851,852; SCS 804; SMK151,153,154,155-1,156,157,SMK181,182,SMV 104T,107,111,168;SMS316;SMC108;S MV104,Servomed N4, Servocard, Serie Cardioserv,Vicom SM,SMU 511,611,612,630 Serie SMK-200, SMK210,211,230,231,240,250,SMP 300,310,Pepita 201,202,202 R.,203,Pepita 203, R204,204R,205, Pepita 205R		S0006NE-L	Exi-Cable, 8ft, Round 10P>DB9F, Suitable for Nellcor probe	077 		C5
3F	COMEN(科曼) :STAR-1000; LANKE (兰科) :LK-8600A; 3F(上海三埃弗): DIONA PMS-8100; 美国CSI雷达:RADO 4200; MAYTAG(美泰克):MTK-2000		S0021BC-L	Exi-Cable, 8ft, Readel 1bt 6P>DB9F, Suitable for BCI probe	014 		B4